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## ART® Implantate

REF

X0006

LOT

2022-03-23

DE

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EN

### ■ INTENDED USE

The ART® spinal fixation system is intended for the posterior stabilization of the thoracic, lumbar and lumbosacral spine.

### ■ PRODUCT DESCRIPTION

The ART® spinal fixation system was developed with the goal of providing the surgeon with a fixation system to carry out dorsal spinal stabilizations simply, quickly and effectively.

Important results of these developments are:

- Polyaxial pedicle screw with rod connector (screw head)
- Pedicle screw with double thread, increasing screw breakage resistance
- Freely rotatable rod connector
- Top-loading system

- Prevention of rod bending owing to the option of adjusting the rod connector by up to  $\pm 30^\circ$  ( $60^\circ$ )

ART® is made from a titanium alloy (Ti6Al4V) in accordance with the standard ASTM F136. ART® must not be used in direct connection with components from other manufacturers. This does not include the use of bone cement.

### ■ SCOPE

These instructions for use are valid for all implants of the ART® product group.

### ■ APPROVED INSTRUMENTS

Only the instruments described in the ART® surgical technique may be used in direct connection with the implant. The surgical technique is available from NGMedical on request.

### ■ INDICATIONS

The ART® spinal fixation system is intended to be used for the following indications:

- to fix and stabilize vertebral segments,
- to support osseous fusion,
- when treating the following acute and chronic instabilities and deformities of the thoracic, lumbar and sacral spine: degenerative spondylolisthesis with objective evidence of a relevant spinal stenosis, fractures, dislocations, scoliosis, kyphosis, spinal tumors and failed prior fusions (pseudoarthrosis),
- when treating severe spondylolisthesis (grade 3 and 4) of the lumbar and lumbosacral vertebral segments among patients with a mature skeletal system.

Furthermore, ART® can also be used in connection with MOVE®-P. Please refer to the MOVE®-P instructions for use and surgical technique.

### ■ CONTRAINDICATIONS

- Obesity from grade II (BMI value  $\geq 35$ ),
- active systemic infection,
- allergy to the implant material (Ti6Al4V),
- endocrine or metabolic illnesses in the patient's history which influence the bone and mineral metabolism,
- health conditions which are opposed to the potential benefit of a spinal operation,
- the patient has received medication influencing bone and mineral metabolism less than two weeks before the planned date of operation,
- bone anomalies which prevent the secure fixation of screws (e.g. dysplastic pedicles),
- the patient is not willing or able to follow the post-operative instructions,
- pregnancy,
- spinal metastases (exception: additional surgical stabilization measures).

ART® was neither developed for indications other than those stated above, nor is it marketed for these.

### ■ COMPLICATIONS AND POSSIBLE ADVERSE EFFECTS

Before the operation, possible adverse effects and further surgery which may be necessary to correct these effects should be explained to the patient:

- adjacent segment degeneration,
- leakage of cerebrospinal fluid,
- failed fusion,
- foreign body reactions to components (caused by allergies), their debris or their fragments,
- infections, including those of the meninges,
- loosening, detachment or breakage of components,
- nerve damage and/or spinal cord damage including the loss of neurological functions (sensory and/or motor deficits such as hypoesthesia, anesthesia, dysesthesia, hyperesthesia, paresthesia, muscular paresis or paralysis, loss of reflexes, increased or reduced muscular tonus, spasticity, bladder, bowel and sexual dysfunctions),
- necessity of further surgery,
- spinal cord injuries with partial or complete paraplegia,
- death,
- visceral or vascular injuries,
- wound healing disorders.

### ■ WARNINGS

The safety and effectiveness of the ART® implants has been proven only for the spinal disorders stated in the indications section. As with all operations, the use of this system involves the potential risk of death. Other potential risks which may require additional surgery are as follows:

- allergic reaction or sensitivity to the titanium alloy (Ti6Al4V),
- deformation and/or breakage of system components,
- infections,
- loosening of components,
- pseudoarthrosis,
- damage to the spinal cord and/or nerves.

By implantation of metals and alloys into the human body, these are exposed to aggressive chemical milieus of salts, acids and alkalis, which might lead to corrosion. If different metals are positioned in close vicinity, galvanic effects may accelerate the corrosion process. In particular, ART® must not be used with components (exception: bone substitute material or fillers) or instruments from other manufacturers.

When removing the implant from the packaging, it must be ensured that the nominal size of the implant matches the size stated on the packaging label.

A completely satisfactory outcome cannot be achieved in every operation. This also applies to spinal operations, where the outcome is determined by many factors.

The osteosynthesis with ART® must normally be supported by an arthrodesis of the affected vertebrae.

### ■ PRECAUTIONS

The surgeon must be completely familiar with ART®, the required instruments and the general surgical technique. The surgeon must use the instruments which are specifically approved by NGMedical for use with ART®. For maximum stability, the surgeon must select the correct size for the pedicle. The steps of the operation can also be read in the surgical technique, which can be requested if needed.

Please document the used implants for each patient with REF, LOT and UDI so that the level of traceability required by law is guaranteed. The implants are intended for single use. An implant must not be reused once explanted.

The ART® implants have not been tested with regard to their safety and compatibility with MRI. They have also not been tested for heating or migration in an MRI environment.

### ■ PACKAGING

ART® implants are generally supplied in double sterile packaging. Only use the implants if the packaging is intact. The storage conditions must guarantee the intact nature of the packaging and therefore also the implants. Defective implants or other components must not be used. Each packaging has a documentation label for you to use and an implantation ID. The packaging and labels contain all the required information. There, you can find the product number REF, LOT, description, method of sterilization, quantity, expiration date and manufacturer as well as the UDI code. All ART® implants are sterilized using gamma radiation.

### ■ COMPLAINTS

Every specialist (e.g. doctor who uses the product) who has a complaint or who is dissatisfied with the quality, labeling, reliability, safety, effectiveness and/or performance of the product should inform NGMedical or, where applicable, the distributor. In the case of a severe incident or the risk of a severe incident which may lead or has led to a severe health impairment or the death of the patient, a user or third-party, NGMedical (or the distributor) and the competent authority must be informed immediately either verbally or in writing. All complaints must be made, stating the product name, REF, LOT, UDI and description as well as the date of the incident. The reporting party must state their name and contact information and describe the incident in as much detail as possible. The legal reporting obligations must also be observed!

### ■ DISPOSAL

For disposal of NGMedical implants, please comply with the specific hospital regulations or return the processed, sterilized and labeled implants to NGMedical.

### ■ EXPECTED LIFETIME


The expected lifetime of ART is designed for 2 years after implantation without fusion. After successful fusion, the expected lifetime corresponds to the remaining lifetime of the patient. In direct combination with MOVE-P, the expected lifetime is 10 years. After exceeding the expected lifetime, a revision may become necessary.

### ■ FURTHER INFORMATION

For further information, please contact:

NGMedical GmbH  
Talsperrenblick 5  
66620 Nonnweiler  
GERMANY  
mail@ngmedical.de  
www.ngmedical.de  
+49 (0) 6873 99997-0

The Summary of Safety and Clinical Performance (SSCP) can be found under the following link: [www.ngmedical.de/sscp/art](http://www.ngmedical.de/sscp/art).

Access to the public website of the Eudamed database: <https://ec.europa.eu/tools/eudamed>.

|                |   |                  |                              |
|----------------|---|------------------|------------------------------|
|                | Consult instructions for use  | <b>STERILE R</b> | Sterilized using irradiation |
|                | Use by  |                  | Manufacturer                 |
|                | Do not re-use   |                  | Storage temperature range    |
|                | Do not resterilize  |                  | Keep dry                     |
| <b>REF</b>     | Catalogue or model number   | <b>LOT</b>       | Batch code                   |
|                | Do not use if package is damaged  |                  | Date of manufacture          |
| <b>MD</b>      | Medical Device  | <b>UDI</b>       | Unique Device Identification |
| <b>Rx ONLY</b> | The medical device according to 21 CFR 801.109(b)(1) may only be prescribed by a physician in accordance with the law (USA) |                  |                              |

© The NGMedical brands and registered trademarks are protected by national and international rights in Germany, the European Union and other regions and countries around the world.

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