GEBRAUCHSANWEISUNG / INSTRUCTIONS FOR USE / MANUEL D'INSTRUCTIONS / ISTRUZIONI PER L'USO /
INSTRUÇÕES DE UTILIZAÇÃO / INSTRUCCIONES DE USO / 使用说明书 / ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ
**BEE® Implantate / BEE® Implants / Implants BEE® / Impianti BEE® / Implante BEE® / Implantes BEE® / BEE® 植入体 /
Εμφυτεύματα BEE®**

DE

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REFX0002LOT2025-01-29

EN | NOT FOR US

■ INTENDED USE
BEE Cages were developed as an intercorporeal implant for the anterior cervical spondylolysis (C3-C7).
■ PRODUCT DESCRIPTION
BEE Cages were especially adapted to the local anatomy to secure the surgical outcome at the best. The caudal side is flat, the cranial side is domed and the implant is formed conically from anterior to posterior. In the lateral view, the implant has a slightly lordotic form. The cage disposes of 4 pins on both bone facing surfaces, which serve the primary fixation. BEE implants are made of titanium alloy (Ti6Al4V) according to ASTM F 3001. BEE Cages are produced with selective laser melting (SLM) technique. It is not allowed to use BEE implants in contact with components of other manufacturers with the exception of bone graft substitutes or fillers.
■ SCOPE

This Instructions for Use is valid for all cervical BEE implants made of titanium alloy.
■ APPROVED INSTRUMENTS
Only the instruments described in the BEE surgical technique (REF X1003 (German), X1004 (English), X1007 (US)), may be used in direct connection with the implant.
■ INDICATIONS
BEE is intended to be used for the following indications:

<ul style="list-style-type: none">• Discopathy,• Disc hernia with radiculopathy and / or myelopathy,• Stenosis of the foramen and spinal canal with or without cervical radiculopathy/myelopathy.
■ CONTRAINDICATIONS

It is not allowed to implant BEE in patients with the following health status:
<ul style="list-style-type: none">• Have a bone mineral density with T-score ≤ -2.5 as determined by spine DXA if male ≥ 60 years of age or female ≥ 50 years of age (Exception: additional operational stabilization measures),• Have an active systemic infection or infection at the surgical site,• Have sustained an osteoporotic fracture / osteopenia of the spine, hip or wrist (Exception: additional operational stabilization measures),• Have received medications (e.g., methotrexate, alendronate) that interfere with bone and mineral metabolism within 2 weeks of the planned date of the index surgery,• Have any medical or surgical condition precluding the potential benefit of spinal surgery,• Have a history of endocrine or metabolic disorders (e.g., Paget's disease) known to affect bone and mineral metabolism,• Have spinal metastases (Exception: additional operational stabilization measures),• Have a known allergy to titanium alloy,• Are pregnant,• Have axial neck pain as the solitary symptom,• Have advanced cervical anatomical deformity (e.g., ankylosing spondylitis, scoliosis) at the operative site,• Are unwilling to follow the postoperative instructions.
■ COMPLICATIONS AND POSSIBLE UNDESIRE OCCURENCES

Prior to the surgery, patients should be explained the possible undesired occurrences and the possible further surgeries which might be necessary for the correction of these:
<ul style="list-style-type: none">• Breakage / failure of system components (loosening, collapse or breaking of components,• Dislocation / migration of the implant,• Possible wound healing disorder and/or wound complications,• Infection,• Neural damages and/or spinal cord damages including loss of neurological functions (sensory and/ or motor deficits such as hypoesthesia, anaesthesia, dyesthesia, hyperaesthesia, paraesthesia, muscle paresis or paralysis, reflex deficiency, increased or decreased muscle tone, spasticity),• Spinal cord injuries with partial or complete transectional lesion,• Leakage of the cerebrospinal fluid,• Spinal fluid infection (myelo-meningitis)• Trachea and/or esophagus injury,• Vocal cord paresis,• Hoarseness,• Foreign body reaction to components or fragments/debris,• Visceral or vascular injuries (Carotid arteries, jugular veins),• Modification of the spinal curvature with kyphotic or scoliotic deformity, loss of correction,• Loss of bone density by stress-shielding,• Need of further surgeries,• Failed Fusion with pseudarthrosis,• Subsidence of the implant into the endplates with loss of intervertebral space height,• Accelerated degeneration of adjacent segments (Adjacent Segment Disease),• Permanent disability,• Osteolysis, bone loss or fracture, bone resorption above or below the level of the surgical segment,• Dysphagia,• Death.
■ WARNINGS

The safety and efficacy of the BEE implants has only been proven for the spinal diseases listed in the paragraph „indications“. There is the potential risk of death when using this system. Further potential risks, which might lead to further surgeries, are:
<ul style="list-style-type: none">• Breakage/Failure of system components,• Loss of Anchorage,• Pseudarthrosis,• Allergical reaction or sensitivity on titanium alloy,• Nerve damage and/or spinal cord damage,• Vascular lesion.

By implantation of metals and alloys into the human body, these are exposed to aggressive chemical milieus of electrolytes, acids and alkalis, which might lead to corrosion. The contact proximity of different metals might accelerate the corrosion process due to galvanic corrosion/reaction. From metallurgical, mechanical and functional reasons, the combination of implant components from different manufacturers is not recommended. Particularly, it is not allowed to use the BEE Cage with components (except for bone graft substitutes or fillers) or surgical techniques of other manufacturers. Generally, implanted materials may cause immune responses or foreign-body reactions with chronic inflammatory changes. It must be ensured that the size of the cage itself matches with the size indicated on the packaging label, when taking the device of its packaging box.

■ PRECAUTIONS
The surgeon must be perfectly familiar with BEE, the necessary instruments as well as the general surgical technique. The surgeon must exclusively use the instruments especially approved by NGMedical for BEE. The surgeon always has to ensure that the optimum implant size is used and that there is an anatomically correct implant placement. Trial implants enable a simple and safe size selection. Over-distraction of the target segment is to be avoided. The surgeon has to ensure the correct implant placement in the cranial, caudal, lateral, anterior and posterior position. The general surgical procedure is described in detail the surgical technique brochure provided by NGMedical (REF X1003 (German), X1004(English), X1007 (US)). Implanted devices are to be documented patient-related with REF and LOT, so that the legally required traceability is guaranteed. The implants are intended for single use. An explanted implant must never be reused.
■ MRI COMPATIBILITY
BEE implants are conditionally MRI compatible. Patients with these implants can be scanned under 1.5 Tesla and 3.0 Tesla static magnetic field strength.
■ PACKAGING
BEE implants are packed double sterile per item. Please only use implants if the packaging box is undamaged. Storage conditions must guarantee the integrity of the packaging box and the implants respectively. Defective implants or further components are not be used. Each box contains documentation labels for your disposal and an implant card. The box and the labels contain all information required. They indicate the REF, LOT, description, sterilization method, quantity, expiration date, UDI-Code and manufacturer. All BEE implants are sterilized by gamma radiation.
■ COMPLAINTS
Any healthcare professional (e. g., a surgeon using the product), who has a complaint or who is dissatisfied with the quality, identification, reliability, safety, efficacy and/or performance of the BEE products, should notify NGMedical or, where applicable, the distributor. In the event of a severe incident or risk of a severe incident, liable to result in or to

have resulted in the death or severe deterioration in the state of health of a patient or user or other persons, NGMedical (or the distributor) and the competent authority has to be informed as soon as possible in writing or verbal. All complaints must include the product name, REF, LOT and the description as well as all data of the incident. The reporting person should state his name, contact data and a detailed description of the incident. The legal reporting obligations are also to be observed!

■ DISPOSAL
For disposal of NGMedical implants, please comply with the specific hospital regulations or return the processed, sterilized and labeled implants to NGMedical.
■ EXPECTED LIFETIME
The expected lifetime of BEE is designed for 2.5 years after implantation under non-fused conditions. After successful fusion, the expected lifetime corresponds to the remaining lifetime of the patient. If the fusion is not successful, a revision may become necessary.
■ FURTHER INFORMATION
Please see chapter "FURTHER INFORMATION" in the next section.

EN | FOR US ONLY

■ CAUTION
Federal law (USA) restricts these devices to sale by or on the order of a physician. These devices should be implanted only by a physician who is fully trained with the devices, intended use, instrumentation and with knowledge of the surgical techniques required. Contact your representative for surgical technique.
■ INTENDED USE
BEE Cages were developed as an intercorporeal implant for the anterior cervical spondylolysis.
■ PRODUCT DESCRIPTION
BEE Cages were especially adapted to the local anatomy to secure the surgical outcome at the best. The caudal side is flat, the cranial side is domed and the implant is formed conically from anterior to posterior. In the lateral view, the implant has a slightly lordotic form. The cage disposes of 4 pins on both bone facing surfaces. BEE implants are made of titanium alloy (Ti6Al4V) according to ASTM-F 3001. It is not allowed to use BEE implants in contact with components of other manufacturers.
■ SCOPE
This Instructions for Use is valid for all cervical BEE implants made of titanium alloy.
■ APPROVED INSTRUMENTS
Only the instruments described in the BEE surgical technique (REF X1007 (US)), may be used in direct connection with the implant.
■ INDICATIONS
BEE Cages are intended for intervertebral body fusion devices in skeletally mature patients for the treatment of cervical disc degeneration and/or cervical spinal instability as confirmed by imaging studies (radiographs, CT, MRI) that results in radiculopathy, myelopathy and/or pain at one or more contiguous levels from C2-T1. These patients should have had at least six weeks of nonoperative treatment. BEE Cages are to be used with autogenous and/or allogeneic bone graft comprised of cancellous and/or corticocancellous bone graft to facilitate fusion and in combination with supplemental fixation indicated for cervical fusion procedures.
■ CONTRAINDICATIONS

It is not allowed to implant BEE in patients with the following health status:
<ul style="list-style-type: none">• Have a bone mineral density with T-score ≤ -2.5 as determined by spine DXA if male ≥ 60 years of age or female ≥ 50 years of age (Exception: additional operational stabilization measures),• Have an active systemic infection or infection at the surgical site,• Have sustained an osteoporotic fracture / osteopenia of the spine, hip or wrist (Exception: additional operational stabilization measures),• Have received medications (e.g., methotrexate, alendronate) that interfere with bone and mineral metabolism within 2 weeks of the planned date of the index surgery,• Have any medical or surgical condition precluding the potential benefit of spinal surgery,• Have a history of endocrine or metabolic disorders (e.g., Paget's disease) known to affect bone and mineral metabolism,• Have spinal metastases (Exception: additional operational stabilization measures),• Have a known allergy to titanium alloy,• Be pregnant,• Have axial neck pain as the solitary symptom,• Have advanced cervical anatomical deformity (e.g., ankylosing spondylitis, scoliosis) at the operative site,• Patient is unwilling to follow the postoperative instructions.• Active infections or high risk of infection,• Signs of local inflammation.• Fever or leukocytosis.• Morbid obesity.• Mental illness.• Any patient having inadequate tissue coverage over the operative site.• Use of this type of surgical implant surgery in children or pediatric patients presents particular risks because of bone growth or physical movement. Subsequent re-intervention may be required.• Infants with a known hereditary or acquired bone friability or calcification problem, or those with a very short life expectancy, should not be considered for this type of surgery.• Substance abuse or senility that precludes the patient from following post operative precautions to prevent implant failure.• Any patient in which implant utilization would interfere with anatomical structures or expected physiological performance.
■ COMPLICATIONS AND POSSIBLE UNDESIRE OCCURENCES

Prior to the surgery, patients should be explained the possible undesired occurrences and the possible further surgeries which might be necessary for the correction of these:
<ul style="list-style-type: none">• Breakage / failure of system components (loosening, collapse or breaking of components,• The risk of migration of the implant components is higher without the use of supplemental fixation,• Possible degeneration of the skin and/or wound complications,• Infection,• Nerval damages and/or spinal cord damages including loss of neurological functions (sense or mobility paralysis), dysaesthesia, hyperaesthesia, paraesthesia, radiculopathy, reflex deficiency,• Spinal cord injuries, outflow of the cerebrospinal fluid,• trachea and/or esophagus injury,• Foreign body reaction to components or fragments/debris,• Visceral or vascular injuries,• Modification of the spinal curvature. Loss of correction, height and/or repositioning,• Loss of bone density by stress-shielding,• Need of further surgeries,• Failed Fusion,• Subsidence of the implant into the endplates,• Adjacent Segment Disease,• Osteolysis, bone loss or fracture, bone resorption above or below the level of the surgical segment,• Dysphagia,• Death.

GEBRAUCHSANWEISUNG / INSTRUCTIONS FOR USE / MANUEL D'INSTRUCTIONS / ISTRUZIONI PER L'USO /
INSTRUÇÕES DE UTILIZAÇÃO / INSTRUCCIONES DE USO / 使用说明书 / ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ
BEE® Implantate / BEE® Implants / Implants BEE® / Impianti BEE® / Implante BEE® / Implantes BEE® / BEE® 植入体 / Εμφυτεύματα BEE®

■ WARNINGS

The safety and efficacy of the BEE implants has only been proven for the spinal diseases listed in the paragraph „indications“. There is the potential risk of death when using this system. Further potential risks, which might lead to a further surgery, are:

- Breakage/Failure of system components,
- Loss of Anchorage,
- Pseudarthrosis,
- Allergic reaction or sensitivity to titanium alloy,
- Nerve damage and/or spinal cord damage,
- Vascular lesion.

By the implantation of metals and alloys into the human body, these are exposed to aggressive chemical milieus of salts, acids and alkalis, which might lead to corrosion. The contact proximity of different metals might accelerate the corrosion process due to galvanic corrosion/reaction. From metallurgical, mechanical and functional reasons, the combination of implant components from different manufacturers is not recommended. It is especially not recommended to use components made from dissimilar metals in direct combination. Particularly, it is not allowed to use the BEE Cage with components or surgical techniques of other manufacturers. It must be ensured that the size of the cage itself matches with the size indicated on the packaging label, when taking the cage out of its packaging box.

■ PRECAUTIONS

The surgeon must be perfectly familiar with BEE, the necessary instruments as well as the general surgical technique. The surgeon must exclusively use the instruments especially approved by NGMedical for BEE. The surgeon always has to ensure that the optimum implant size is used and that there is an anatomically correct implant placement. Trial implants enable a simple and safe size destination. Over distraction is to be avoided. The surgeon must pay attention that there is a correct implant placement in the cranial and caudal as well as in the anterior and posterior position. The general surgical procedure can be learned out of the surgical technique (REF X1007 (US)). This document can be obtained from NGMedical or the representative.

Single use only. Implants are provided as single use only implants and are not to be reused or implanted regardless of an apparent undamaged condition.

Please document the used implants patient-related with REF and LOT, so that the legally required traceability is guaranteed.

Based on the fatigue testing results, the physician/surgeon should consider the levels of implantation, patient weight, patient activity level, other patient conditions, etc. which may impact on the performance of the system.

■ MRI COMPATIBILITY

BEE has not been evaluated for safety and compatibility in the MR environment. It has not been tested for heating migration, or image artifact in the MR environment. The safety of the BEE system in the MR environment is unknown. Scanning a patient who has this device may result in patient injury.

■ PACKAGING

BEE implants are packed double sterile per item. Please only use implants if the packaging box is undamaged. Storage conditions must guarantee the integrity of the packaging box and the implants respectively. Defective implants or further components are not to be used. Each box contains documentation labels for your disposal and an implant card. The box and the labels contain all information required. They indicate the REF, LOT, description, sterilization method, quantity, expiration date, UDI-Code and manufacturer. All BEE implants are sterilized by gamma radiation.

■ COMPLAINTS

Any healthcare professional (e. g., a surgeon using the product), who has a complaint or who is dissatisfied with the quality, identification, reliability, safety, efficacy and/or performance of the BEE products, should notify NGMedical or, where applicable, the distributor. In the event of a severe incident or risk of a severe incident, liable to result in or to have resulted in the death or severe deterioration in the state of health of a patient or user or other persons, NGMedical (or the distributor) has to be informed as soon as possible in writing or verbal. All complaints must include the product name, REF, LOT and the description as well as all data of the incident. The reporting person should state his name, contact data and a detailed description of the incident. The legal reporting obligations are also to be observed!

■ DISPOSAL


For disposal of NGMedical implants, please comply with the specific hospital regulations or return the processed, sterilized and labeled implants to NGMedical.

■ FURTHER INFORMATION

For further information please contact

NGMedical GmbH
Trierer Straße 65
66620 Nonnweiler
GERMANY
mail@ngmedical.de
www.ngmedical.de
+49 (0) 6873 99997-0

The Summary of Safety and Clinical Performance (SSCP) can be found under the following link:
www.ngmedical.de/sscp/bee.

Access to the public website of the Eudamed database: <https://ec.europa.eu/tools/eudamed>.

Symbol	Standard	Ref no	Title	Meaning
	ISO 15223-1- Symbols to be used with information to be supplied by the manufacturer - Part 1: General requirements	5.4.3	Consult instructions for use	Indicates the need for the user to consult the instructions for use.
		5.1.4	Use by	Indicates the date after which the medical device is not to be used.
		5.4.2	Do not re-use	Indicates a medical device that is intended for one use, or for use on a single patient during a single procedure.
		5.2.6	Do not re-sterilize	Indicates a medical device that is not to be re-sterilized.
		5.1.6	Catalogue or model number	Indicates the manufacturer's catalogue number so that the medical device can be identified.
		5.2.8	Do not use if package is damaged	Indicates a medical device that should not be used if the package has been damaged or opened.
		5.2.4	Sterilized using irradiation	Indicates a medical device that has been sterilized using irradiation.
		5.1.1	Manufacturer	Indicates the medical device manufacturer.
		5.3.7	Storage temperature range	Indicates the temperature limits to which the medical device can be safely exposed.
		5.3.4	Keep dry	Indicates a medical device that needs to be protected from moisture.

		5.1.5	Batch code	Indicates the manufacturer's batch code so that the batch or lot can be identified.
		5.1.3	Date of manufacture	Indicates the date of manufacture the medical device
		5.7.7	Medical Device	Indicates the item is a medical device
		5.7.10	Unique Device Identification	Indicates a carrier that contains Unique Device Identifier information
	Double sterile barrier system			
	European conformity mark with Notified Body identification number			
	The medical device according to 21 CFR 801.109(b)(1) may only be prescribed by a physician in accordance with the law (USA)			

® The NGMedical brands and registered trademarks are protected by national and international rights in Germany, the European Union and other regions and countries worldwide.

GEBRAUCHSANWEISUNG / INSTRUCTIONS FOR USE / MANUEL D'INSTRUCTIONS / ISTRUZIONI PER L'USO /
INSTRUÇÕES DE UTILIZAÇÃO / INSTRUCCIONES DE USO / 使用说明书 / ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ
BEE® Implantate / BEE® Implants / Implants BEE® / Impianti BEE® / Implante BEE® / Implantes BEE® / BEE® 植入体 / Εμφυτεύματα BEE®

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© Die NGMedical Marken und eingetragenen Handelsmarken sind durch nationale und internationale Rechte in Deutschland, der Europäischen Union und anderen Regionen und Ländern weltweit geschützt.