GEBRAUCHSANWEISUNG / INSTRUCTIONS FOR USE / MANUEL D'INSTRUCTIONS / ISTRUZIONI PER L'USO /
INSTRUÇÕES DE UTILIZAÇÃO / INSTRUCCIONES DE USO / 使用说明书 / ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ / NÁVOD K POUŽITÍ /
NÁVOD NA POUŽITIE / HASZNÁLATI UTASÍTÁS / INSTRUCȚIUNI DE UTILIZARE
Instrumente / Instruments / Instruments / Strumenti / Instrumentos / Instrumentos / 医疗器械 / Εργαλεία / Nástroje /
Nástroje / Hangszerek / Instrumente

DE

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EN | NOT FOR US

1. INTENDED PURPOSE

The instruments of class Ir are reusable surgical instruments with the purpose of enabling the implantation of an implant from NGMedical. The surgical technique of each implant describes the instruments to be used in each case.

NGMedical's trials are reusable, surgically invasive products for temporary use. They are used for sizing as well as for determining the final position of the implant.

The surgical technique of each implant describes the instruments to be used. The instruments may only be used in direct connection with implants from NGMedical and as described in the respective surgical techniques.

2. SCOPE

These instructions for use apply to all NGMedical instruments. The instruments of NGMedical are intended for surgical use by appropriately trained surgical staff.

3. IMPLANTS

All surgical techniques are available from NGMedical on request. The article numbers of the implants can be found in the corresponding surgical techniques.

4. INDICATIONS/CONTRAINDICATIONS

The NGMedical instruments are only used in direct connection with the associated implant. Therefore, the indications/contraindications of the respective implant apply. The indications/contraindications of the associated implant are as follows:

MOVE®-C

INDICATIONS

MOVE®-C is intended to be used for the following indications:

- Discopathy,
- Disc hernia with radiculopathy,
- Stenosis of the foramen and spinal canal.

CONTRAINDICATIONS

MOVE®-C should not be implanted in patients with the following conditions:

- Bone mineral density with T-score ≤ -1.5 as determined by spine DXA if male ≥ 60 years of age or female ≥ 50 years of age,
- Active systemic infection or infection at the operative site,
- Sustained osteoporotic fracture of the spine, hip or wrist,
- Medications (e.g., methotrexate, alendronate) that interfere with bone and mineral metabolism within 2 weeks of the planned date of the index surgery,
- Any medical or surgical condition precluding the potential benefit of spinal surgery,
- History of endocrine or metabolic disorders (e.g., Paget's disease) known to affect bone and mineral metabolism,
- Spinal metastases,
- Known allergy to titanium (Ti6Al4V_ELI) or polycarbonate-urethane,
- Pregnancy,
- Axial neck pain as the solitary symptom,
- Severe cervical myelopathy as evidenced by any sign of gait disturbance, unilateral or bilateral leg weakness, and/or uncontrollable bowel/bladder symptoms related to cervical spine disease,
- Treatment (e.g., posterior element decompression) that destabilizes the spine,
- Advanced cervical anatomical deformity (e.g., ankylosing spondylitis, scoliosis) at the operative site,
- Advanced degenerative changes (e.g., spondylosis) at the index vertebral level as evidenced by:
 - Bridging osteophytes,
 - Average ROM $< 4^\circ$,
 - Disc height $< 25\%$ of the height of the inferior vertebral disc space, as measured in a lateral radiograph in neutral position,
 - Subluxation > 3 mm,
 - Kyphotic deformity at $> 20^\circ$ on neutral radiographs.

During the postoperative rehabilitation, there is a certain risk for patients who are partly or not fully in the position to follow the necessary precautions. Factors like the degree of physical activity and the following of the instructions regarding the postoperative treatment affect the stress, the implant is exposed to. There is a higher pseudarthrosis-risk for smokers. These patients are to be informed on the risk and to be requested to stop smoking prior to and immediately after surgery.

BEE®

INDICATIONS

BEE® is intended to be used for the following indications:

- Discopathy,
- Disc hernia with radiculopathy and / or myelopathy,
- Stenosis of the foramen and spinal canal with or without cervical radiculopathy/myelopathy.

CONTRAINDICATIONS

It is not allowed to implant BEE® in patients with the following health status:

- Bone mineral density with T-score ≤ -2.5 as determined by spine DXA if male ≥ 60 years of age or female ≥ 50 years of age (Exception: additional operational stabilization measures),

- Active systemic infection or infection at the surgical site,
- Sustained osteoporotic fracture / osteopenia of the spine, hip or wrist (Exception: additional operational stabilization measures),
- Medications (e.g., methotrexate, alendronate) that interfere with bone and mineral metabolism within 2 weeks of the planned date of the index surgery,
- Any medical or surgical condition precluding the potential benefit of spinal surgery,
- History of endocrine or metabolic disorders (e.g., Paget's disease) known to affect bone and mineral metabolism,
- Spinal metastases (Exception: additional operational stabilization measures),
- Known allergy to titanium alloy (Ti6Al4V_ELI),
- Pregnancy,
- Axial neck pain as the solitary symptom,
- Advanced cervical anatomical deformity (e.g., ankylosing spondylitis, scoliosis) at the operative site,
- The patient is not willing or not able to follow the postoperative instructions.

BEE® HA

INDICATIONS

BEE® HA is intended to be used for the following indications:

- Discopathy,
- Disc hernia with radiculopathy and / or myelopathy,
- Stenosis of the foramen and spinal canal with or without cervical radiculopathy/ myelopathy.

CONTRAINDICATIONS

It is not allowed to implant BEE® HA in patients with the following health status:

- Bone mineral density with T-score ≤ -2.5 as determined by spine DXA if male ≥ 60 years of age or female ≥ 50 years of age (Exception: additional operational stabilization measures),
- Active systemic infection or infection at the surgical site,
- Sustained osteoporotic fracture/ osteopenia of the spine, hip or wrist (Exception: additional operational stabilization measures),
- Medications (e.g., methotrexate, alendronate) that interfere with bone and mineral metabolism within 2 weeks of the planned date of the index surgery,
- Any medical or surgical condition precluding the potential benefit of spinal surgery,
- History of endocrine or metabolic disorders (e.g., Paget's disease) known to affect bone and mineral metabolism,
- Spinal metastases (Exception: additional operational stabilization measures),
- Known allergy to polyetheretherketone (PEEK), hydroxyapatite (HA) or titanium alloy (Ti6Al4V_ELI),
- Pregnancy,
- Axial neck pain as the solitary symptom,
- Advanced cervical anatomical deformity (e.g., ankylosing spondylitis, scoliosis) at the operative site,
- The patient not willing or able to follow the postoperative instructions.

ART®

INDICATIONS

The ART® spinal fixation system is intended to be used for the following indications:

- To fix and stabilize vertebral segments,
- To support osseous fusion,
- When treating the following acute and chronic instabilities and deformities of the thoracic, lumbar and sacral spine: degenerative spondylolisthesis with objective evidence of a relevant spinal stenosis, fractures, dislocations, scoliosis, kyphosis, spinal tumors and failed prior fusions (pseudoarthrosis),
- When treating severe spondylolisthesis (grade 3 and 4) of the lumbar and lumbosacral vertebral segments among patients with a mature skeletal system.

Furthermore, ART® can also be used in connection with MOVE®-P. Please refer to the MOVE®-P instructions for use and surgical technique.

CONTRAINDICATIONS

- Obesity from grade II (BMI value ≥ 35),
- Active systemic infection,
- Allergy to the implant material (Ti6Al4V_ELI),
- Endocrine or metabolic illnesses in the patient's history which influence the bone and mineral metabolism,
- Health conditions which are opposed to the potential benefit of a spinal operation,
- The patient has received medication influencing bone and mineral metabolism less than two weeks before the planned date of operation,
- Bone anomalies which prevent the secure fixation of screws (e.g. dysplastic pedicles),
- The patient is not willing or able to follow the post-operative instructions,
- Pregnancy,
- Spinal metastases (exception: additional surgical stabilization measures).

ART® was neither developed for indications other than those stated above, nor is it marketed for these.

BEE® PLIF and NGM WAVE® PLIF

INDICATIONS

BEE® PLIF and NGM WAVE® PLIF are intended for the following indications:

- Degenerative disorders of the spine including degenerative disc disease, vertebral canal and foraminal stenosis,
- Lumbar spine instability in one or more segments,
- Pseudarthrosis,
- Spondylolisthesis, including degenerative spondylolisthesis,
- Spondylodiscitis.

CONTRAINDICATIONS

BEE® PLIF and NGM WAVE® PLIF must not be implanted in patients with the following health conditions:

- Obesity from grade II (BMI value ≥ 35),
- Active systemic infection,
- Allergy to the implant material (see the imprint on the packaging),
- Endocrine or metabolic illness in the patient's history which influences the bone and mineral metabolism,
- Health conditions which are opposed to the potential benefit of a spinal surgery,
- The patient has received medication influencing bone and mineral metabolism within two weeks prior to the planned date of operation,
- Osteoporosis,
- The patient is not willing or able to follow the post-operative instructions,
- Pregnancy,
- Vertebral fracture(s),
- Spinal metastases (exception: additional surgical stabilization measures).

All instructions for use are available from NGMedical on request.

5. POSSIBLE COMPLICATIONS

Possible complications are:

- Loss of function of the instrument,
- Laser marking becomes illegible,
- Corrosion,
- Damaged surfaces,
- Split-offs,
- Contamination,
- Fissures and sharp edges,
- Deformation,
- Breakage/failure of components (loosening, disintegration or breakage of components).

Sort out the instruments indicating one or several complications. All instruments must be checked for their functionality according to the surgical technique before the start of an operation. Please also note the document "Reprocessing Instructions Instruments" (REF X0200).

6. TRAYS

Each instrument system of NGMedical has its own tray with lid, in which the instruments can be inserted.

GEBRAUCHSANWEISUNG / INSTRUCTIONS FOR USE / MANUEL D'INSTRUCTIONS / ISTRUZIONI PER L'USO /
INSTRUÇÕES DE UTILIZAÇÃO / INSTRUCCIONES DE USO / 使用说明书 / ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ / NÁVOD K POUŽITÍ /
NÁVOD NA POUŽITIE / HASZNÁLATI UTASÍTÁS / INSTRUCȚIUNI DE UTILIZARE
Instrumente / Instruments / Instruments / Strumenti / Instrumentos / Instrumentos / 医疗器械 / Εργαλεία / Nástroje /
Nástroje / Hangszerek / Instrumente

REF X0008 LOT 2025-11-03

7. REPROCESSING

All NGMedical instruments are placed on the market non-sterile and must be subjected to reprocessing and sterilisation before each use in accordance with the document "Reprocessing Instructions Instruments" (REF X0200). The Reprocessing Instructions Instruments are available from NGMedical on request.

8. COMPLAINTS

Any healthcare professional (e.g., a surgeon using the product), who has a complaint or who is dissatisfied with the quality, identification, reliability, safety, efficacy and/or performance of the products, should notify NGMedical or, where applicable, the distributor. In the event of a severe incident or risk of a severe incident, liable to result in or to have resulted in the death or severe deterioration in the state of health of a patient or user or other persons, NGMedical (or the distributor) and the competent authority has to be informed as soon as possible in writing or verbal. All complaints must include the product name, REF, SN and the description as well as all data of the incident. The reporting person should state his name, contact data and a detailed description of the incident. The legal reporting obligations are also to be observed!

9. DISPOSAL

NGMedical's instruments should be disposed of in accordance with the specific hospital disposal guidelines or the reprocessed, sterilised and accordingly labelled instruments can be sent to NGMedical.

10. FURTHER INFORMATION

For further information, please contact:

NGMedical GmbH
Trierer Straße 65
66620 Nonnweiler
GERMANY
mail@ngmedical.de
www.ngmedical.de
+49 (0) 6873 99997-0

CE/CE 0482

	Consult Instructions for Use		Non-Sterile
	Keep Dry		Manufacturer
	Reference-No.		Serial Number
	UDI (Unique Device Identification)		Medical Device
	Do not Use if Package is Damaged		Date of Manufacture
	European Conformity Mark		European Conformity Mark with Notified Body Identification Number

© The NGMedical brands and registered trademarks are protected by national and international rights in Germany, the European Union and other regions and countries worldwide.

EN | FOR US ONLY

1. INTENDED PURPOSE

The instruments of class I_r are reusable surgical instruments with the purpose of enabling the implantation of an implant from NGMedical. The surgical technique of each implant describes the instruments to be used in each case.

NGMedical's trials are reusable, surgically invasive products for temporary use. They are used for sizing as well as for determining the final position of the implant.

The surgical technique of each implant describes the instruments to be used. The instruments may only be used in direct connection with implants from NGMedical and as described in the respective surgical techniques.

2. SCOPE

These instructions for use apply to all NGMedical instruments. The instruments of NGMedical are intended for surgical use by appropriately trained surgical staff.

3. IMPLANTS

All surgical techniques are available from NGMedical on request. The article numbers of the implants can be found in the corresponding surgical techniques.

4. INDICATIONS/CONTRAINDICATIONS

The NGMedical instruments are only used in direct connection with the associated implant. Therefore, the indications/contraindications of the respective implant apply. The indications/contraindications of the associated implant can be viewed in the respective instructions for use.

MOVE®-C

INDICATIONS

MOVE®-C is intended to be used for the following indications:

- Discopathy,
- Disc hernia with radiculopathy,
- Stenosis of the foramen and spinal canal.

CONTRAINDICATIONS

MOVE®-C should not be implanted in patients with the following conditions:

- Bone mineral density with T-score ≤ -1.5 as determined by spine DXA if male ≥ 60 years of age or female ≥ 50 years of age,
- Active systemic infection or infection at the operative site,
- Sustained osteoporotic fracture of the spine, hip or wrist,
- Medications (e.g., methotrexate, alendronate) that interfere with bone and mineral metabolism within 2 weeks of the planned date of the index surgery,
- Any medical or surgical condition precluding the potential benefit of spinal surgery,
- History of endocrine or metabolic disorders (e.g., Paget's disease) known to affect bone and mineral metabolism,
- Spinal metastases,
- Known allergy to titanium (Ti6Al4V_ELI) or polycarbonate-urethane,
- Pregnancy,
- Axial neck pain as the solitary symptom,
- Severe cervical myelopathy as evidenced by any sign of gait disturbance, unilateral or bilateral leg weakness, and/or uncontrollable bowel/bladder symptoms related to cervical spine disease,
- Treatment (e.g., posterior element decompression) that destabilizes the spine,
- Advanced cervical anatomical deformity (e.g., ankylosing spondylitis, scoliosis) at the operative site,
- Advanced degenerative changes (e.g., spondylosis) at the index vertebral level as evidenced by:
 - Bridging osteophytes,
 - Average ROM $< 4^\circ$,
 - Disc height $< 25\%$ of the height of the inferior vertebral disc space, as measured in a lateral radiograph in neutral position,
 - Subluxation > 3 mm,

- Kyphotic deformity at $> 20^\circ$ on neutral radiographs.

During the postoperative rehabilitation, there is a certain risk for patients who are partly or not fully in the position to follow the necessary precautions. Factors like the degree of physical activity and the following of the instructions regarding the postoperative treatment affect the stress, the implant is exposed to. There is a higher pseudarthrosis-risk for smokers. These patients are to be informed on the risk and to be requested to stop smoking prior to and immediately after surgery.

BEE®

INDICATIONS

BEE® Cages are intended for intervertebral body fusion devices in skeletally mature patients for the treatment of cervical disc degeneration and/or cervical spinal instability as confirmed by imaging studies (radiographs, CT, MRI) that results in radiculopathy, myelopathy and/or pain at one or more contiguous levels from C2-T1. These patients should have had at least six weeks of nonoperative treatment. BEE® Cages are to be used with autogenous and/or allogeneic bone graft comprised of cancellous and/or corticocancellous bone graft to facilitate fusion and in combination with supplemental fixation indicated for cervical fusion procedures.

CONTRAINDICATIONS

It is not allowed to implant BEE® in patients with the following health status:

- Bone mineral density with T-score ≤ -2.5 as determined by spine DXA if male ≥ 60 years of age or female ≥ 50 years of age (Exception: additional operational stabilization measures),
- Active systemic infection or infection at the surgical site,
- Sustained osteoporotic fracture / osteopenia of the spine, hip or wrist (Exception: additional operational stabilization measures),
- Medications (e.g., methotrexate, alendronate) that interfere with bone and mineral metabolism within 2 weeks of the planned date of the index surgery,
- Any medical or surgical condition precluding the potential benefit of spinal surgery,
- History of endocrine or metabolic disorders (e.g., Paget's disease) known to affect bone and mineral metabolism,
- Spinal metastases (Exception: additional operational stabilization measures),
- Known allergy to titanium alloy (Ti6Al4V_ELI),
- Pregnancy,
- Axial neck pain as the solitary symptom,
- Advanced cervical anatomical deformity (e.g., ankylosing spondylitis, scoliosis) at the operative site,
- The patient is not willing or able to follow the postoperative instructions,
- Active infections or high risk of infection,
- Signs of local inflammation,
- Fever or leukocytosis,
- Morbid obesity,
- Mental illness,
- Any patient having inadequate tissue coverage over the operative site,
- Use of this type of surgical implant surgery in children or pediatric patients presents particular risks because of bone growth or physical movement. Subsequent re-intervention may be required,
- Infants with a known hereditary or acquired bone friability or calcification problem, or those with a very short life expectancy, should not be considered for this type of surgery,
- Substance abuse or senility that precludes the patient from following post operative precautions to prevent implant failure,
- Any patient in which implant utilization would interfere with anatomical structures or expected physiological performance.

BEE® HA

INDICATIONS

BEE® HA Cages are intended for intervertebral body fusion devices in skeletally mature patients for the treatment of cervical disc degeneration and/or cervical spinal instability as confirmed by imaging studies (radiographs, CT, MRI) that results in radiculopathy, myelopathy and/or pain at one or more contiguous levels from C2-T1. These patients should have had at least six weeks of nonoperative treatment. BEE® HA Cages are to be used with autogenous and/or allogeneic bone graft comprised of cancellous and/or corticocancellous bone graft to facilitate fusion and in combination with supplemental fixation indicated for cervical fusion procedures.

CONTRAINDICATIONS

It is not allowed to implant BEE® HA in patients with the following health status:

- Bone mineral density with T-score ≤ -2.5 as determined by spine DXA if male ≥ 60 years of age or female ≥ 50 years of age (Exception: additional operational stabilization measures),
- Active systemic infection or infection at the surgical site,
- Sustained osteoporotic fracture/ osteopenia of the spine, hip or wrist (Exception: additional operational stabilization measures),
- Medications (e.g., methotrexate, alendronate) that interfere with bone and mineral metabolism within 2 weeks of the planned date of the index surgery,
- Any medical or surgical condition precluding the potential benefit of spinal surgery,
- History of endocrine or metabolic disorders (e.g., Paget's disease) known to affect bone and mineral metabolism,
- Spinal metastases (Exception: additional operational stabilization measures),
- Known allergy to polyetheretherketone (PEEK), hydroxyapatite (HA) or titanium alloy (Ti6Al4V_ELI),
- Pregnancy,
- Axial neck pain as the solitary symptom,
- Advanced cervical anatomical deformity (e.g., ankylosing spondylitis, scoliosis) at the operative site,
- The patient is not willing or able to follow the postoperative instructions,
- Active infections or high risk of infection,
- Signs of local inflammation,
- Fever or leukocytosis,
- Morbid obesity,
- Mental illness,
- Any patient having inadequate tissue coverage over the operative site,
- Use of this type of surgical implant surgery in children or pediatric patients presents particular risks because of bone growth or physical movement. Subsequent re-intervention may be required,
- Infants with a known hereditary or acquired bone friability or calcification problem, or those with a very short life expectancy, should not be considered for this type of surgery,
- Substance abuse or senility that precludes the patient from following post operative precautions to prevent implant failure,
- Any patient in which implant utilization would interfere with anatomical structures or expected physiological performance.

ART®

INDICATIONS

The ART® spinal fixation system is intended to be used for the following indications:

- To fix and stabilize vertebral segments,
 - To support osseous fusion,
 - When treating the following acute and chronic instabilities and deformities of the thoracic, lumbar and sacral spine: degenerative spondylolisthesis with objective evidence of a relevant spinal stenosis, fractures, dislocations, scoliosis, kyphosis, spinal tumors and failed prior fusions (pseudoarthrosis),
 - When treating severe spondylolisthesis (grade 3 and 4) of the lumbar and lumbosacral vertebral segments among patients with a mature skeletal system.
- Furthermore, ART® can also be used in connection with MOVE®-P. Please refer to the MOVE®-P instructions for use and surgical technique.

CONTRAINDICATIONS

- Obesity from grade II (BMI value ≥ 35),
- Active systemic infection,
- Allergy to the implant material (Ti6Al4V_ELI),
- Endocrine or metabolic illnesses in the patient's history which influence the bone and mineral metabolism,
- Health conditions which are opposed to the potential benefit of a spinal operation,
- The patient has received medication influencing bone and mineral metabolism less than two weeks before the planned date of operation,
- Bone anomalies which prevent the secure fixation of screws (e.g. dysplastic pedicles),
- The patient is not willing or able to follow the post-operative instructions,

GEBRAUCHSANWEISUNG / INSTRUCTIONS FOR USE / MANUEL D'INSTRUCTIONS / ISTRUZIONI PER L'USO /
INSTRUÇÕES DE UTILIZAÇÃO / INSTRUCCIONES DE USO / 使用说明书 / ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ / NÁVOD K POUŽITÍ /
NÁVOD NA POUŽITIE / HASZNÁLATI UTASÍTÁS / INSTRUCȚIUNI DE UTILIZARE
Instrumente / Instruments / Instruments / Strumenti / Instrumentos / Instrumentos / 医疗器械 / Εργαλεία / Nástroje /
Nástroje / Hangszerek / Instrumente

REF X0008 LOT 2025-11-03

- Pregnancy,
 - Spinal metastases (exception: additional surgical stabilization measures).
- ART® was neither developed for indications other than those stated above, nor is it marketed for these.

BEE® PLIF and NGM WAVE® PLIF

INDICATIONS

BEE® PLIF and NGM WAVE® PLIF are intended for the following indications:

- Degenerative disorders of the spine including degenerative disc disease, vertebral canal and foraminal stenosis,
- Lumbar spine instability in one or more segments,
- Pseudarthrosis,
- Spondylolisthesis, including degenerative spondylolisthesis,
- Spondylodiscitis.

CONTRAINDICATIONS

BEE® PLIF and NGM WAVE® PLIF must not be implanted in patients with the following health conditions:

- Obesity from grade II (BMI value ≥ 35),
- Active systemic infection,
- Allergy to the implant material (see the imprint on the packaging),
- Endocrine or metabolic illness in the patient's history which influences the bone and mineral metabolism,
- Health conditions which are opposed to the potential benefit of a spinal surgery,
- The patient has received medication influencing bone and mineral metabolism within two weeks prior to the planned date of operation,
- Osteoporosis,
- The patient is not willing or able to follow the post-operative instructions,
- Pregnancy,
- Vertebral fracture(s),
- Spinal metastases (exception: additional surgical stabilization measures).

All instructions for use are available from NGMedical on request.

5. POSSIBLE COMPLICATIONS

Possible complications are:

- Loss of function of the instrument,
- Laser marking becomes illegible,
- Corrosion,
- Damaged surfaces,
- Split-offs,
- Contamination,
- Fissures and sharp edges,
- Deformation,
- Breakage/failure of components (loosening, disintegration, or breakage of components).

Sort out the instruments indicating one or several complications. All instruments must be checked for their functionality according to the surgical technique before the start of an operation. Please also note the document "Reprocessing Instructions Instruments" (REF X0200).

6. TRAYS

Each instrument system of NGMedical has its own tray with lid, in which the instruments can be inserted.

7. REPROCESSING

All NGMedical instruments are placed on the market non-sterile and must be subjected to reprocessing and sterilization before each use in accordance with the document "Reprocessing Instructions Instruments" (REF X0200). The Reprocessing Instructions Instruments are available from NGMedical on request.

8. DISPOSAL

NGMedical's instruments should be disposed of in accordance with the specific hospital disposal guidelines or the reprocessed, sterilized and accordingly labelled instruments can be sent to NGMedical.

9. FURTHER INFORMATION

For further information, please contact:

NGMedical GmbH
Trierer Straße 65
66620 Nonnweiler
GERMANY
mail@ngmedical.de
www.ngmedical.de
+49 (0) 6873 99997-0

CE/CE 0482

	Consult Instructions for Use		Non-Sterile
	Keep Dry		Manufacturer
	Reference-No.		Serial Number
	UDI (Unique Device Identification)		Medical Device
	Do not use if Package is Damaged		Date of Manufacture
	European Conformity Mark		European Conformity Mark with Notified Body Identification Number

© The NGMedical brands and registered trademarks are protected by national and international rights in Germany, the European Union and other regions and countries worldwide.

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