GEBRAUCHSANWEISUNG / INSTRUCTIONS FOR USE / MANUEL D'INSTRUCTIONS / ISTRUZIONI PER L'USO / INSTRUÇÕES DE UTILIZAÇÃO / INSTRUCCIONES DE USO / 使用说明书 / ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ / NÁVOD K POUŽITÍ / NÁVOD NA POUŽITIE / HASZNÁLATI UTASÍTÁS / INSTRUȚIUNI DE UTILIZARE

MOVE®-C Implantate / Implants / Implants MOVE®-C / Impianti MOVE®-C / Implantes MOVE®-C / Implantes MOVE®-C / MOVE®-C 植入体 / Εμφυτεύματα MOVE®-C / Implantáty MOVE®-C / Implantát MOVE®-C / MOVE®-C implantátumok / Implanturi MOVE®-C

DE

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EN | NOT FOR US

■ INTENDED PURPOSE

The MOVE®-C cervical disc prosthesis is a prosthesis for the dynamic stabilization of one or two cervical segments (C3-C7) after anterior discectomy has been effected by taking into consideration that the movability of the segment to be treated with the prosthesis is to be maintained at its best for patients with a fully-grown skeletal system.

■ PRODUCT DESCRIPTION

The MOVE®-C implant is delivered mounted on an implant holder and consists of a cranial titanium plate (coated with titanium niobium nitride (TiNbN) on the inner surface) and a caudal titanium plate (Ti6Al4V\_EL1) with moulded polymer core out of PCU (Polycarbonate – urethane). The plates possess pins on the bone facing surfaces, which serve the primary fixation. The implant holder consists of PEEK (Polyetheretherketone).  
MOVE®-C implants have been especially adapted to the anatomy in order to secure the surgical result at the best. The caudal endplate is flat, the cranial endplate is domed and the implant is formed conically from anterior to posterior. In the lateral view, the implant has a slightly lordotic form.  
It is not allowed to use MOVE®-C implants in direct combination with components of other manufacturers.

■ SCOPE

This Instructions for Use is valid for all MOVE®-C implants.

■ INTENDED USERS

These instructions for use alone do not provide sufficient information for immediate use of the system. Training in the use of the products by an experienced and trained surgeon and/or by NGMedical (or the distributor) is strongly recommended.  
The procedure must be performed in accordance with the instructions for use and the surgical technique. The surgeon is responsible for ensuring that the operation is performed correctly. It is strongly recommended that the procedure is only performed by surgeons who are qualified and experienced in spinal surgery, are aware of the general risks associated with spinal surgery and are familiar with the product-specific surgical techniques.

■ APPROVED INSTRUMENTS

Only the instruments described in the MOVE®-C surgical techniques may be used in direct connection with the implant (Please see the respective IFU REF X0008). The surgical technique is available from NGMedical on request.

■ EXPECTED CLINICAL BENEFIT

The MOVE®-C cervical disc prosthesis is intended for use in accordance with its intended purpose and instructions for use for the dynamic stabilization of adjacent segments of the cervical spine. The expected clinical benefit of surgery with MOVE®-C is the preservation of physiological mobility of the treated segment in all 6 degrees of freedom. It is intended to relieve pain, restore function to the affected segment, and reduce the risk of adjacent segment degeneration. The clinical benefit is supported by preclinical data and study data.  
The Summary of Safety and Clinical Performance (SSCP) can be found under the following link:  
[www.ngmedical.de/sscp/move\\_c](http://www.ngmedical.de/sscp/move_c).

Access to the public website of the Eudamed database: <https://ec.europa.eu/tools/eudamed>.

■ INDICATIONS

MOVE®-C is intended to be used for the following indications:

- Discopathy,
- Disc hernia with radiculopathy,
- Stenosis of the foramen and spinal canal.

■ CONTRAINDICATIONS

MOVE®-C should not be implanted in patients with the following conditions:

- Bone mineral density with T-score ≤ -1.5 as determined by spine DXA if male ≥ 60 years of age or female ≥ 50 years of age,
- Active systemic infection or infection at the operative site,
- Sustained osteoporotic fracture of the spine, hip or wrist,
- Medications that interfere with bone and mineral metabolism within 2 weeks of the planned date of the index surgery,
- Any medical or surgical condition precluding the potential benefit of spinal surgery,
- History of endocrine or metabolic disorders known to affect bone and mineral metabolism,
- Spinal metastases,
- Known allergy to titanium or polycarbonate-urethane,
- Pregnancy,
- Axial neck pain as the solitary symptom,
- Severe cervical myelopathy as evidenced by any sign of gait disturbance, unilateral or bilateral leg weakness, and/or uncontrollable bowel/bladder symptoms related to cervical spine disease,
- Required treatment (e.g., posterior element decompression) that destabilizes the spine,
- Advanced cervical anatomical deformity (e.g., ankylosing spondylitis, scoliosis) at the operative site,
- Advanced degenerative changes (e.g., spondylosis) at the index vertebral level as evidenced by:
  - Bridging osteophytes,
  - Average ROM < 4°,
  - Disc height < 25% of the height of the inferior vertebral disc space, as measured in a lateral radiograph in neutral position,
  - Subluxation > 3 mm,
  - Kyphotic deformity at > 20° on neutral radiographs.

During the postoperative rehabilitation, there is a certain risk for patients who are partly or not fully in the position to follow the necessary precautions. Factors like the degree of physical activity and the following of the instructions regarding the postoperative treatment affect the stress, the implant is exposed to. There is a higher pseudarthrosis-risk for smokers. These patients are to be informed on the risk and to be requested to stop smoking prior to and immediately after surgery.

■ COMPLICATIONS AND POSSIBLE UNDESIRED OCCURENCES

Prior to the surgery, patients should be explained the possible undesired occurrences and the possible further surgeries which might be necessary for the correction of these:

- Loosening, collapse or breaking of components,
- Migration of the implant components,
- Sensitive reaction of the soft tissue on the implant material,
- Possible degeneration of the skin and/or wound complications,
- Pseudarthrosis or reluctant ingrowth or bad ingrowth,
- Infection,
- Nerval damages and/or spinal cord damages including loss of neurological functions (sense or mobility paralysis), dysaesthesia, hyperaesthesia, paraesthesia, radiculopathy, reflex deficiency,
- Spinal cord injuries, outflow of the cerebrospinal fluid,
- Trachea and/or esophagus injury,
- Vocal cord paresis,
- Hoarseness,
- Vertebral fracture,
- Foreign body reaction (allergically based) to components or fragments/debris,
- Loss of anchorage, (e.g.: through strong extension before the MOVE®-C implant has attached to the end plates)
- Visceral or vascular injuries,
- Modification of the spinal curvature. Loss of correction, height and/or repositioning,
- Loss of bone density by stress-shielding,
- Osteolysis, bone loss or fracture, bone resorption above or below the level of the surgical segment,
- Impairment of activities,
- Lack of effective means for the treatment of the symptoms, which the surgery was intended to fix,
- Need of further surgeries,
- Death.

■ WARNINGS

The safety and efficacy of the MOVE®-C implants has only been proven for the spinal diseases listed in the paragraph „Indications“. There is the potential risk of death when using this system. Further potential risks, which might lead to a further surgery are:

- Breakage/Failure of system components,
- Loss of Anchorage, (e.g.: through strong extension before the MOVE®-C implant has attached to the end plates)
- Pseudarthrosis,
- Allergical reaction or sensitivity on metals or PCU,
- Vertebral fracture,
- Nerve damage and/or spinal cord damage,
- Vascular lesion.

By the implantation of metals and alloys into the human body, these are exposed to aggressive chemical milieus of salts, acids and alkalis, which might lead to corrosion. The contact proximity of different metals might accelerate the corrosion process due to galvanic corrosion/reaction. From metallurgical, mechanical and functional reasons, the combination of implant components from different manufacturers is not recommended. Particularly, it is not allowed to use the MOVE®-C disc prosthesis with components or surgical techniques of other manufacturers. It must be ensured that the size of the prosthesis itself matches with the size indicated on the packaging label, when taking the prosthesis out of its packaging box.

■ PRECAUTIONS

The surgeon must be perfectly familiar with MOVE®-C, the required instruments as well as the general surgical technique. The surgeon must exclusively use the instruments especially approved by NGMedical for MOVE®-C. The surgeon always has to ensure that the optimum implant size is used and that there is an anatomically correct implant placement. Trial implants enable a simple and safe size selection. Over distraction is to be avoided. The surgeon has to ensure that there is a correct implant placement in the cranial and caudal as well as in the anterior and posterior position. The general surgical procedure is described in detail in the surgical technique provided by NGMedical (REF X1001).

Please document the used implants patient-related with REF, LOT and UDI, so that the legally required traceability is guaranteed. The implants are intended for single use. **An explanted implant must never be reused.**

■ MRI COMPATIBILITY

MOVE®-C implants have been tested for the effects of magnetic resonance imaging (MRI) and can be classified as conditionally MRI-compatible. Patients with these implants can be scanned at static magnetic field strengths of 1.5 Tesla and 3.0 Tesla. The scan duration should not exceed 15 minutes.

■ PACKAGING

MOVE®-C implants are generally supplied in double sterile packaging. Please only use implants if the packaging box is undamaged. Storage conditions must guarantee the integrity of the packaging box and the implants respectively. Defective implants or other components must not be used. Each packaging box contains documentation labels for free use and an implant card. The packaging box and the labels contain all information required. They indicate REF, LOT, description, sterilization method, quantity, expiration date, manufacturing date UDI-Code and manufacturer. All MOVE®-C implants are sterilized using gamma radiation or ethylene oxide gas. The type of sterilization is indicated on each packaging box.

■ COMPLAINTS

Any healthcare professional (e.g., a surgeon using the product), who has a complaint or who is dissatisfied with the quality, identification, reliability, safety, efficacy and/or performance of the MOVE®-C products, should notify NGMedical or, where applicable, the distributor. In the event of a severe incident or risk of a severe incident, liable to result in or to have resulted in the death or severe deterioration in the state of health of a patient or user or other persons, NGMedical (or the distributor) and the competent authority must be informed as soon as possible in writing or verbal. All complaints must include the product name, REF, LOT, UDI and the description as well as all data of the incident. The reporting person should state his name, contact data and a detailed description of the incident. The legal reporting obligations are also to be observed!

■ DISPOSAL


For disposal of NGMedical implants, please comply with the specific hospital regulations or return the processed, sterilized and labeled implants to NGMedical.

■ EXPECTED LIFETIME

The expected lifetime of MOVE®-C is set at 15 years after implantation. After exceeding the expected lifetime, a revision may become necessary.

■ FURTHER INFORMATION









For further information please contact

	NGMedical GmbH Trierer Straße 65 66620 Nonnweiler GERMANY mail@ngmedical.de <a href="http://www.ngmedical.de">www.ngmedical.de</a> +49 (0) 6873 99997-0			
	 			
Symbol	Standard	Ref no	Title	Meaning
	ISO 15223-1- Symbols to be used with information to be supplied by the	5.4.3	Consult instructions for use	Indicates the need for the user to consult the instructions for use.

GEBRAUCHSANWEISUNG / INSTRUCTIONS FOR USE / MANUEL D'INSTRUCTIONS / ISTRUZIONI PER L'USO / INSTRUÇÕES DE UTILIZAÇÃO / INSTRUCCIONES DE USO / 使用说明书 / ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ / NÁVOD K POUŽITÍ / NÁVOD NA POUŽITIE / HASZNÁLATI UTASÍTÁS / INSTRUȚIUNI DE UTILIZARE

MOVE®-C Implantate / Implants / Implants MOVE®-C / Impianti MOVE®-C / Implantes MOVE®-C / Implantes MOVE®-C / MOVE®-C 植入体 / Εμφυτεύματα MOVE®-C / Implantáty MOVE®-C / Implantát MOVE®-C / MOVE®-C implantátumok / Implanturi MOVE®-C

REF X0001 LOT 2025-08-06

	manufacturer - Part 1: General requirements	5.1.4	Use by	Indicates the date after which the medical device is not to be used.
		5.4.2	Do not re-use	Indicates a medical device that is intended for one use, or for use on a single patient during a single procedure.
		5.2.6	Do not resterilize	Indicates a medical device that is not to be resterilized.
		5.1.6	Catalogue or model number	Indicates the manufacturer's catalogue number so that the medical device can be identified.
		5.2.8	Do not use if package is damaged	Indicates a medical device that should not be used if the package has been damaged or opened.
		5.2.4	Sterilized using irradiation	Indicates a medical device that has been sterilized using irradiation.
		5.2.3	Sterilized with ethylene oxide	Indicates a medical device that was sterilized with ethylene oxide
		5.1.1	Manufacturer	Indicates the medical device manufacturer.
		5.3.7	Storage temperature range	Indicates the temperature limits to which the medical device can be safely exposed.
		5.3.4	Keep dry	Indicates a medical device that needs to be protected from moisture.
		5.1.5	Batch code	Indicates the manufacturer's batch code so that the batch or lot can be identified.
		5.1.3	Date of manufacture	Indicates the date of manufacture the medical device
		5.7.7	Medical Device	Indicates the item is a medical device
		5.7.10	Unique Device Identification	Indicates a carrier that contains Unique Device Identifier information
	Double sterile barrier system			
	European conformity mark with Notified Body identification number			
	The medical device according to 21 CFR 801.109(b)(1) may only be prescribed by a physician in accordance with the law (USA)			
	Conditionally MRI compatible			

© The NGMedical brands and registered trademarks are protected by national and international rights in Germany, the European Union and other regions and countries worldwide.

EN | FOR US ONLY

CAUTION

Federal law (USA) restricts these devices to sale by or on the order of a physician. These devices should be implanted only by a physician who is fully trained with the devices, intended use, instrumentation and with knowledge of the surgical techniques required. Contact your representative for surgical technique.

INTENDED PURPOSE

The MOVE®-C cervical disc prosthesis is a prosthesis for the dynamic stabilization of one or two cervical segments (C3-C7) after anterior discectomy has been effected by taking into consideration that the movability of the segment to be treated with the prosthesis is to be maintained at its best for patients with a fully-grown skeletal system.

PRODUCT DESCRIPTION

The MOVE®-C implant is delivered mounted on an implant holder and consists of a cranial titanium plate (coated with titanium niobium nitride (TiNbN) on the inner surface) and a caudal titanium plate (Ti6Al4V\_ELJ) with moulded polymer core out of PCU (Polycarbonate – urethane). The plates possess pins on the bone facing surfaces, which serve the primary fixation. The implant holder consists of PEEK (Polyetheretherketone).

MOVE®-C implants have been especially adapted to the anatomy in order to secure the surgical result at the best. The caudal endplate is flat, the cranial endplate is domed and the implant is formed conically from anterior to posterior. In the lateral view, the implant has a slightly lordotic form.

It is not allowed to use MOVE®-C implants in direct combination with components of other manufacturers.

SCOPE

This Instructions for Use is valid for all MOVE®-C implants.

APPROVED INSTRUMENTS

All instruments of the product group MOVE®-C (Please see the respective IFU (REF X0008)).

INDICATIONS

MOVE®-C is intended to be used for the following indications:

- Discopathy,
- Disc hernia with radiculopathy,
- Stenosis of the foramen and spinal canal.

CONTRAINDICATIONS

MOVE®-C should not be implanted in patients with the following conditions:

- Bone mineral density with T-score  $\leq -1.5$  as determined by spine DXA if male  $\geq 60$  years of age or female  $\geq 50$  years of age,
- Active systemic infection or infection at the operative site,
- Sustained osteoporotic fracture of the spine, hip or wrist,
- Medications that interfere with bone and mineral metabolism within 2 weeks of the planned date of the index surgery,
- Any medical or surgical condition precluding the potential benefit of spinal surgery,
- History of endocrine or metabolic disorders known to affect bone and mineral metabolism,
- Spinal metastases,
- Known allergy to titanium or polycarbonate-urethane,
- Pregnancy,
- Axial neck pain as the solitary symptom,
- Severe cervical myelopathy as evidenced by any sign of gait disturbance, unilateral or bilateral leg weakness, and/or uncontrollable bowel/bladder symptoms related to cervical spine disease,
- Required treatment (e.g., posterior element decompression) that destabilizes the spine,
- Advanced cervical anatomical deformity (e.g., ankylosing spondylitis, scoliosis) at the operative site,
- Advanced degenerative changes (e.g., spondylosis) at the index vertebral level as evidenced by:
  - Bridging osteophytes,
  - Average ROM  $< 4^\circ$ ,
  - Disc height  $< 25\%$  of the height of the inferior vertebral disc space, as measured in a lateral radiograph in neutral position,
  - Subluxation  $> 3$  mm,
  - Kyphotic deformity at  $> 20^\circ$  on neutral radiographs.

During the postoperative rehabilitation, there is a certain risk for patients who are partly or not fully in the position to follow the necessary precautions. Factors like the degree of physical activity and the following of the instructions regarding the postoperative treatment affect the stress, the implant is exposed to. There is a higher pseudarthrosis-risk for smokers. These patients are to be informed on the risk and to be requested to stop smoking prior to and immediately after surgery.

COMPLICATIONS AND POSSIBLE UNDESIED OCCURENCES

Prior to the surgery, patients should be explained the possible undesired occurrences and the possible further surgeries which might be necessary for the correction of these:

- Loosening, collapse or breaking of components,
- Migration of the implant components,
- Sensitive reaction of the soft tissue on the implant material,
- Possible degeneration of the skin and/or wound complications,
- Pseudarthrosis or reluctant ingrowth or bad ingrowth,
- Infection,
- Nerval damages and/or spinal cord damages including loss of neurological functions (sense or mobility paralysis), dysaesthesia, hyperaesthesia, paraesthesia, radiculopathy, reflex deficiency,
- Spinal cord injuries, outflow of the cerebrospinal fluid,
- Trachea and/or esophagus injury,
- Vocal cord paresis,
- Hoarseness,
- Vertebral fracture,
- Foreign body reaction (allergically based) to components or fragments/debris,
- Loss of anchorage, (e.g.: through strong extension before the MOVE®-C implant has attached to the end plates)
- Visceral or vascular injuries,
- Modification of the spinal curvature. Loss of correction, height and/or repositioning,
- Loss of bone density by stress-shielding,
- Osteolysis, Loss of bone or fracture, bone resorption above or below the height of the surgical segment,
- Impairment of activities,
- Lack of effective means for the treatment of the symptoms, which the surgery was intended to fix,
- Need of further surgeries,
- Death.

WARNINGS

The safety and efficacy of the MOVE®-C implants has only been proven for the spinal diseases listed in the paragraph „Indications“. There is the potential risk of death when using this system. Further potential risks, which might lead to a further surgery are:

- Breakage/Failure of system components,
- Loss of Anchorage, (e.g.: through strong extension before the MOVE®-C implant has attached to the end plates)
- Pseudarthrosis,
- Allergical reaction or sensitivity on metals or PCU,
- Vertebral fracture,
- Nerve damage and/or spinal cord damage,
- Vascular lesion.

By the implantation of metals and alloys into the human body, these are exposed to aggressive chemical milieus of salts, acids and alkalis, which might lead to corrosion. The contact proximity of different metals might accelerate the corrosion process due to galvanic corrosion/reaction. From metallurgical, mechanical and functional reasons, the combination of implant components from different manufacturers is not recommended. Particularly, it is not allowed to use the MOVE®-C disc prosthesis with components or surgical techniques of other manufacturers. It must be ensured that the size of the prosthesis itself matches with the size indicated on the packaging label, when taking the prosthesis out of its packaging box.

PRECAUTIONS

The surgeon must be perfectly familiar with MOVE®-C, the required instruments as well as the general surgical technique. The surgeon must exclusively use the instruments especially approved by NGMedical for MOVE®-C. The surgeon always has to ensure that the optimum implant size is used and that there is an anatomically correct implant placement. Trial implants enable a simple and safe size selection. Over distraction is to be avoided. The surgeon has to ensure that there is a correct implant placement in the cranial and caudal as well as in the anterior and posterior position. The general surgical procedure is described in detail in the surgical technique provided by NGMedical (REF X1001).

Please document the used implants patient-related with REF, LOT and UDI, so that the legally required traceability is guaranteed. The implants are intended for single use. **An explanted implant must never be reused.**

MRI COMPATIBILITY

MOVE® implants have been tested for the effects of magnetic resonance imaging (MRI) and can be classified as conditionally MRI-compatible. Patients with these implants can be scanned at static magnetic field strengths of 1.5 Tesla and 3.0 Tesla. The scan duration should not exceed 15 minutes.

PACKAGING

MOVE®-C implants are generally supplied in double sterile packaging. Please only use implants if the packaging box is undamaged. Storage conditions must guarantee the integrity of the packaging box and the implants respectively. Defective implants or other components must not be used. Each packaging box contains documentation labels for free use and an implant card. The packaging box and the labels contain all information required. They indicate REF, LOT, description, sterilization method, quantity, expiration date, manufacturing date UDI code and manufacturer. All MOVE®-C-implants are sterilized by gamma radiation.

COMPLAINTS

Any healthcare professional (e.g., a surgeon using the product), who has a complaint or who is dissatisfied with the quality, identification, reliability, safety, efficacy and/or performance of the MOVE®-C products, should notify

GEBRAUCHSANWEISUNG / INSTRUCTIONS FOR USE / MANUEL D'INSTRUCTIONS / ISTRUZIONI PER L'USO / INSTRUÇÕES DE UTILIZAÇÃO / INSTRUCCIONES DE USO / 使用说明书 / ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ / NÁVOD K POUŽITÍ / NÁVOD NA POUŽITIE / HASZNÁLATI UTASÍTÁS / INSTRUȚIUNI DE UTILIZARE

MOVE®-C Implantate / Implants / Implants MOVE®-C / Impianti MOVE®-C / Implantes MOVE®-C / Implants MOVE®-C / MOVE®-C 植入体 / Εμφυτεύματα MOVE®-C / Implantáty MOVE®-C / Implantát MOVE®-C / MOVE®-C implantátumok / Implanturi MOVE®-C

NGMedical or, where applicable, the distributor. In the event of a severe incident or risk of a severe incident, liable to result in or to have resulted in the death or severe deterioration in the state of health of a patient or user or other persons, NGMedical (or the distributor) and the competent authority must be informed as soon as possible in writing or verbal. All complaints must include the product name, REF, LOT, UDI and the description as well as all data of the incident. The reporting person should state his name, contact data and a detailed description of the incident. The legal reporting obligations are also to be observed!

DISPOSAL

For disposal of NGMedical implants, please comply with the specific hospital regulations or return the processed, sterilized and labeled implants to NGMedical.

EXPECTED LIFETIME

The expected lifetime of MOVE®-C is set at 15 years after implantation. After exceeding the expected lifetime, a revision may become necessary.

FURTHER INFORMATION

For further information please contact

NGMedical GmbH  
Trierer Straße 65  
66620 Nonnweiler  
GERMANY  
mail@ngmedical.de  
www.ngmedical.de  
+49 (0) 6873 99997-0

The Summary of Safety and Clinical Performance (SSCP) can be found under the following link: [www.ngmedical.de/sscp/move\\_c](http://www.ngmedical.de/sscp/move_c).

Access to the public website of the Eudamed database: <https://ec.europa.eu/tools/eudamed>.

© The NGMedical brands and registered trademarks are protected by national and international rights in Germany, the European Union and other regions and countries worldwide.

Symbol	Standard	Ref no	Title	Meaning
	ISO 15223-1- Symbols to be used with information to be supplied by the manufacturer - Part 1: General requirements	5.4.3	Consult instructions for use	Indicates the need for the user to consult the instructions for use.
		5.1.4	Use by	Indicates the date after which the medical device is not to be used.
		5.4.2	Do not re-use	Indicates a medical device that is intended for one use, or for use on a single patient during a single procedure.
		5.2.6	Do not re-sterilize	Indicates a medical device that is not to be re-sterilized.
		5.1.6	Catalogue or model number	Indicates the manufacturer's catalogue number so that the medical device can be identified.
		5.2.8	Do not use if package is damaged	Indicates a medical device that should not be used if the package has been damaged or opened.
		5.2.4	Sterilized using irradiation	Indicates a medical device that has been sterilized using irradiation.
		5.1.1	Manufacturer	Indicates the medical device manufacturer.
		5.3.7	Storage temperature range	Indicates the temperature limits to which the medical device can be safely exposed.
		5.3.4	Keep dry	Indicates a medical device that needs to be protected from moisture.
		5.1.5	Batch code	Indicates the manufacturer's batch code so that the batch or lot can be identified.
		5.1.3	Date of manufacture	Indicates the date of manufacture the medical device
		5.7.7	Medical Device	Indicates the item is a medical device
		5.7.10	Unique Device Identification	Indicates a carrier that contains Unique Device Identifier information
		Double sterile barrier system		
	European conformity mark with Notified Body identification number			
	The medical device according to 21 CFR 801.109(b)(1) may only be prescribed by a physician in accordance with the law (USA)			

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