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## PLIF Implantate

REF X0009 LOT 2022-03-23

DE

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EN

### ■ INTENDED USE

BEE® PLIF cages and NGM WAVE® PLIF cages are implants for intervertebral, posterior, lumbar fusion in the segments L2-S1 among patients with a mature skeletal system.

### ■ PRODUCT DESCRIPTION

BEE® PLIF cages and NGM WAVE® PLIF cages have been specially adapted to suit the local anatomy to ensure the best possible surgical outcome. BEE® PLIF cages are made of titanium alloy in accordance with ASTM F3001 and are produced in an additive manufacturing process. NGM WAVE® PLIF cages are made of PEEK-OPTIMA™ Natural and tantalum (pins) in accordance with ASTM F560. The implants are available in various heights, lengths and lordosis angles. The cranial and caudal side of the implant is rough or fitted with teeth to achieve a good anchoring of the implant.

### ■ SCOPE

These instructions for use are valid for all BEE® PLIF and NGM WAVE® PLIF implants.

### ■ APPROVED INSTRUMENTS

Only the instruments described in the BEE® PLIF and NGM WAVE® PLIF surgical techniques may be used in direct connection with the implant. The surgical technique is available from NGMedical on request.

### ■ INDICATIONS

BEE® PLIF and NGM WAVE® PLIF are intended for the following indications:

- degenerative disorders of the spine including degenerative disc disease, vertebral canal and foraminal stenosis,
- lumbar spine instability in one or several segments,
- pseudoarthrosis,
- spondylolisthesis, including degenerative spondylolisthesis,
- spondylodiscitis.

### ■ CONTRAINDICATIONS

BEE® PLIF and NGM WAVE® PLIF must not be implanted in the following patients or health conditions:

- obesity from grade II (BMI value  $\geq 35$ ),
- active systemic infection,
- allergy to the implant material (see the imprint on the packaging).
- endocrine or metabolic illness in the patient's history which influences the bone and mineral metabolism,
- health conditions which are opposed to the potential benefit of a spinal operation,
- the patient has received medication influencing bone and mineral metabolism less than two weeks before the planned date of operation,
- osteoporosis,
- the patient is unwilling or unable to follow the post-operative instructions,
- pregnancy,
- vertebral fracture(s),
- spinal metastases (exception: additional surgical stabilization measures).

### ■ COMPLICATIONS AND POSSIBLE ADVERSE EFFECTS

Before the operation, the possible adverse effects and possibly required further surgery should be explained to the patient:

- leakage of cerebrospinal fluid,
- spinal curvature with kyphotic or scoliotic deformity, loss of correction,
- (accelerated) adjacent segment degeneration,
- breakage/failure of system components (loosening, degradation or breakage of components),
- permanent disability,
- dislocation/migration of the implant,
- subsidence of the implant into the vertebral body end plates causing reduced height of the intervertebral space,
- failed fusion with pseudoarthrosis,
- foreign body reactions to components, their debris or their fragments,
- infection,
- loss of bone density due to stress shielding,
- possibly required further surgery,
- nerve damage and/or spinal cord damage including the loss of neurological functions (sensory and/or motor deficiencies such as hypoaesthesia, anesthesia, dysesthesia, hyperaesthesia, paresthesia, partial or complete muscular paralysis, lack of reflexes, increased or reduced muscular tonus, spasticity),
- osteoporotic vertebral fracture,
- spinal cord injuries with partial or complete paraplegia,
- death,
- visceral or vascular injuries,
- spinal infections (e.g. spondylitis, discitis, myelo-meningitis)
- wound healing disorders.

### ■ WARNINGS

The safety and effectiveness of the BEE® PLIF and NGM WAVE® PLIF implants have been proven only for the spinal disorders stated in the indications section. As with all operations, the use of this system involves the potential risk of death. Other potential risks which may require further surgery are as follows:

- allergic reaction or sensitivity to the implant material (see the imprint on the packaging).
- breakage/failure of system components,
- vascular injuries,
- pseudoarthrosis,
- damage to the spinal cord and/or nerves,
- loss of anchorage.

By implantation of metals and alloys into the human body, these are exposed to aggressive chemical milieus of salts, acids and alkalis which might lead to corrosion. If different metals are positioned in close vicinity, galvanic effects may accelerate the corrosion process. In particular, BEE® PLIF and NGM WAVE® PLIF must not be used with components (exception: bone substitute material or fillers) or instruments from other manufacturers. In general, implanted material may cause immune reactions or foreign body reactions as well as chronic inflammatory changes.

When removing the implant from the packaging, it must be ensured that the nominal height of the implant matches the size stated on the packaging labeling.

### ■ PRECAUTIONS

The surgeon must be completely familiar with BEE® PLIF and NGM WAVE® PLIF, the required instruments and the general surgical technique. The surgeon must use the instruments which are specifically approved by NGMedical for use with BEE® PLIF and NGM WAVE®. The surgeon must always ensure that the adequate implant size is used and that the implant is fitted into the anatomically correct position. Overdistraction must be avoided. It must be observed that the implant is aligned correctly in the cranial-caudal as well as the ventral-dorsal plane. The steps of the operation can also be read in the implant-specific surgical technique.

Documentation of the used implants is necessary for each patient with REF, LOT and UDI so that the level of traceability required by law is guaranteed. The implants are intended for single use. An implant must not be reused once explanted.

The BEE® PLIF and NGM WAVE® PLIF implants have not been tested with regard to their safety and compatibility with MRI. They have also not been tested for heating or migration in an MRI environment.

### ■ PACKAGING

BEE® PLIF and NGM WAVE® PLIF implants are generally supplied in double sterile packaging. Only use the implants if the packaging is intact. The storage conditions must guarantee the intact nature of the packaging and therefore also the implants. Defective implants must not be used. Each packaging has a documentation label for you to use and an implantation ID. The packaging and labels contain all the required information. There, you can find the product number REF, LOT, description, method of sterilization, quantity, expiration date and manufacturer, as well as scan the UDI code. All BEE® PLIF and NGM WAVE® PLIF implants are sterilized by way of gamma irradiation sterilization.

### ■ COMPLAINTS

Every specialist (e.g. doctor who uses the product) who has a complaint or who is dissatisfied with the quality, labeling, reliability, safety, effectiveness and/or performance of the product should inform NGMedical or, where applicable, the distributor. In the case of a severe incident or the risk of a severe incident which may lead or has led to a severe health impairment or the death of the patient, a user or third-party, NGMedical (or the retailer) and the competent authority must be informed immediately either verbally or in writing. All complaints must be made, stating the product name, REF, LOT, UDI and description as well as the date of the incident. The reporting party must state their name and contact information and describe the incident in as much detail as possible. The legal reporting obligations must also be observed!

### ■ DISPOSAL

For disposal of NGMedical implants, please comply with the specific hospital regulations or return the processed, sterilized and labeled implants to NGMedical.

### ■ EXPECTED LIFETIME


The expected lifetime of BEE® PLIF and NGM WAVE® is designed for 2.5 years after implantation under non-fused conditions. After successful fusion, the expected lifetime corresponds to the remaining lifetime of the patient. If the fusion is not successful, a revision may become necessary.

### ■ FURTHER INFORMATION

For further information, please contact:

NGMedical GmbH  
Talsperrenblick 5  
66620 Nonnweiler  
GERMANY  
mail@ngmedical.de  
www.ngmedical.de  
+49 (0) 6873 99997-0

CE 0482

The Summary of Safety and Clinical Performance (SSCP) can be found under the following link: [www.ngmedical.de/sscp/bee\\_plif](http://www.ngmedical.de/sscp/bee_plif), [www.ngmedical.de/sscp/ngm\\_wave](http://www.ngmedical.de/sscp/ngm_wave)

Access to the public website of the Eudamed database: <https://ec.europa.eu/tools/eudamed>.

	Consult instructions for use	<b>STERILE R</b>	Sterilized using irradiation
	Use by		Manufacturer
	Do not re-use		Storage temperature range
	Do not re-sterilize		Keep dry
<b>REF</b>	Catalogue or model number	<b>LOT</b>	Batch code
	Do not use if package is damaged		Date of manufacture
<b>MD</b>	Medical Device	<b>UDI</b>	Unique Device Identification
<b>Rx ONLY</b>	The medical device according to 21 CFR 801.109(b)(1) may only be prescribed by a physician in accordance with the law (USA)		

© The NGMedical brands and registered trademarks are protected by national and international rights in Germany, the European Union and other regions and countries around the world.

# INSTRUCTIONS FOR USE

## PLIF implants

**REF**      **X0009**      **LOT**      2022-03-23

### EN-FOR US ONLY

#### ■ CAUTION

Federal law (USA) restricts these devices to sale by or on the order of a physician. These devices should be implanted only by a physician who is fully trained with the devices, intended use, instrumentation and with knowledge of the surgical techniques required. Contact your representative for surgical technique.

#### ■ INTENDED USE

BEE® PLIF cages and NGM WAVE® PLIF cages are implants for intervertebral, posterior, lumbar fusion in the segments L2-S1 among patients with a mature skeletal system.

#### ■ PRODUCT DESCRIPTION

BEE® PLIF cages and NGM WAVE® PLIF cages have been specially adapted to suit the local anatomy to ensure the operative outcome to the greatest extent possible. BEE® PLIF cages are made of titanium alloy in accordance with ASTM F3001 and are produced in an additive manufacturing process. NGM WAVE® PLIF cages are made of PEEK-OPTIMA™ Natural and tantalum (pins) in accordance with ASTM F560. The implants are available in various heights, lengths and lordosis angles. The cranial and caudal side of the implant is rough or fitted with teeth to achieve a good degree of anchoring of the implant.

#### ■ SCOPE

These instructions for use are valid for all BEE® PLIF implants and NGM WAVE® PLIF.

#### ■ APPROVED INSTRUMENTS

Only the instruments described in the BEE® PLIF surgical techniques may be used in direct connection with the implant. The removal of the implant is described in the surgical technique, which you can obtain from NGMedical upon request.

#### ■ INDICATIONS

BEE® PLIF cages and NGM WAVE® PLIF cages are indicated for interbody fusion with autogenous bone graft in patients with degenerative disc disease (DDD) at one or two contiguous levels from L2 to S1. These DDD patients may also have up to Grade I spondylolisthesis or retrolisthesis at the involved levels. DDD is defined as discogenic back pain with degeneration of the disc confirmed by history and radiographic studies. These patients should be skeletally mature and have had six months of non-operative treatment. These implants may be implanted via an open or a minimally invasive posterior approach. These implants are to be used with autogenous and/or allogeneic bone graft comprised of cancellous and/or corticocancellous bone graft to facilitate fusion. These devices are intended to be used with supplemental fixation instrumentation, which has been cleared by the FDA for use in the lumbar spine.

#### ■ CONTRAINDICATIONS

BEE® PLIF and NGM WAVE® PLIF cages must not be implanted in the following patients or health situations:

- obesity from grade II (BMI value  $\geq 35$ ),
- active system infection,
- allergy to the implant material (see the imprint on the packaging).
- endocrine or metabolic illness in the patient's history which influences the bone and mineral balance,
- health conditions which are opposed to the potential benefit of a spinal operation,
- the patient took medication which influences the bone and mineral balance less than two weeks before the planned date of operation,
- osteoporosis,
- the patient is not willing or able to follow the post-operative instructions,
- pregnancy,
- vertebral fracture(s),
- spinal metastases (exception: additional surgical stabilization measures).

#### ■ COMPLICATIONS AND POSSIBLE ADVERSE EFFECTS

Before the operation, the possible adverse effects and further surgery which may be necessary to correct the effects should be explained to the patient:

- discharge of cerebrospinal fluid,
- change to the spinal curvature with kyphotic or scoliotic deformity, loss of correction,
- (accelerated) adjacent segment degeneration,
- breakage/failure of system components (loosening, degradation or breakage of components),
- permanent disability,
- dislocation/migration of the implant,
- sinking of the implant into the vertebral body end plates (subsidence) with a reduced height of the intervertebral space,
- failed fusion with pseudoarthrosis,
- foreign body granuloma around components, their debris or their fragments,
- infection,
- loss of bone density due to stress shielding,
- possible further surgery,
- nerve damage and/or spinal cord damage including the loss of neurological functions (sensory and/or motor deficiencies such as hypoaesthesia, anesthesia, dysesthesia, hyperaesthesia, paresthesia, partial or complete muscular paralysis, lack of reflexes, increased or reduced muscular tonus, spasticity),
- osteoporotic vertebral fracture,
- spinal cord injuries with partial or complete paraplegia,
- death,
- visceral or vascular injuries,
- spinal infections (e.g. spondylitis, discitis, myelo-meningitis)
- wound healing disorders.

#### ■ WARNINGS

The safety and effectiveness of the BEE® PLIF and NGM WAVE® PLIF implants have only been proven for the spinal disorders stated in the indications section. As with all operations, the use of this system involves the potential risk of death. Other potential risks which may require further surgery are as follows:

- allergic reaction or sensitivity to the implant material (see the imprint on the packaging).
- breakage/failure of system components,
- vascular injuries,
- pseudoarthrosis,
- damage to the spinal cord and/or nerves,
- loss of anchorage.

When metals and alloys are implanted into the human body, these are exposed to aggressive chemical environments involving salts, acids and bases which can lead to corrosion. If different metals are positioned next to one another in touching proximity, galvanic corrosion may accelerate the corrosion process. In particular, BEE® PLIF and NGM WAVE® PLIF must not be used with components (exception: autogenous and/or allogeneic bone graft comprised of cancellous and/or corticocancellous bone graft to facilitate fusion) or instruments from other manufacturers. In general, implanted material may cause immune reactions or foreign body reactions as well as chronic inflammatory changes. When removing the implant from the packaging, it must be ensured that the nominal height of the implant matches the size stated on the packaging labeling.

#### ■ PRECAUTIONS

The surgeon must be completely familiar with BEE® PLIF and NGM WAVE® PLIF, the required instruments and the general surgical technique. The surgeon must use the instruments which are specifically approved by NGMedical for use with BEE® PLIF and NGM WAVE® PLIF. Based on the dynamic testing results, the physician should consider the levels of implantation, patient weight, patient activity level, other patient conditions, etc., which may impact on the performance of the intervertebral body fusion device. The surgeon must always ensure that the ideal implant size is used and that the implant is fitted into the anatomically correct position. Overdistraction must be avoided. Here, it must be observed that the implant is aligned correctly both in terms of the cranial-caudal line and the ventral-dorsal line. The steps of the operation can also be read in the implant-specific surgical technique.

Please document the used implants for each patient with REF, LOT and UDI so that the level of traceability required by law is guaranteed. The implants are intended for single use. An implant must not be reused once explanted.

#### ■ MRI COMPATIBILITY

BEE® PLIF and NGM WAVE® PLIF have not been evaluated for safety and compatibility in the MR environment. It has not been tested for heating migration, or image artifact in the MR environment. The safety of the BEE® PLIF and NGM WAVE® PLIF systems in the MR environment is unknown. Scanning a patient who has this device may result in patient injury.

#### ■ PACKAGING

BEE® PLIF and NGM WAVE® PLIF implants are generally supplied in double sterile packaging. Only use the implants if the packaging is intact. The storage conditions must guarantee the intact nature of the packaging and therefore also the implants. Defective implants must not be used. Each packaging has a documentation label for you to use and an implantation ID. The packaging and labels contain all the required information. There, you can find the product number REF, LOT, description, method of sterilization, quantity, expiration date and manufacturer, as well as scan the UDI code. All BEE® PLIF and NGM WAVE® PLIF implants are sterilized by way of gamma irradiation sterilization.

#### ■ COMPLAINTS

Every specialist (e.g. doctor who uses the product) who has a complaint or who is dissatisfied with the quality, labeling, reliability, safety, effectiveness and/or performance of the product should inform NGMedical or, where applicable, the

distributor. In the case of a severe incident or the risk of a severe incident which may lead or has led to a severe health impairment or the death of the patient, a user or third-party, NGMedical (or the retailer) must be informed immediately either verbally or in writing. All complaints must be made, stating the product name, REF, LOT, UDI and description as well as the date of the incident. The reporting party must state their name and contact information and describe the incident in as much detail as possible. The legal reporting obligations must also be observed!


#### ■ DISPOSAL

For disposal of NGMedical implants, please comply with the specific hospital regulations or return the processed, sterilized and labeled implants to NGMedical.

#### ■ FURTHER INFORMATION


For further information, please contact:

 NGMedical GmbH  
Talsperrenblick 5  
66620 Nonnweiler  
GERMANY  
mail@ngmedical.de  
www.ngmedical.de  
+49 (0) 6873 99997-0

The Summary of Safety and Clinical Performance (SSCP) can be found under the following link:  
[www.ngmedical.de/sscp/bee\\_plif](http://www.ngmedical.de/sscp/bee_plif), [www.ngmedical.de/sscp/ngm\\_wave](http://www.ngmedical.de/sscp/ngm_wave)

Access to the public website of the Eudamed database: <https://ec.europa.eu/tools/eudamed>.

	Consult instructions for use	<b>STERILE R</b>	Sterilized using irradiation
	Use by		Manufacturer
	Do not re-use		Storage temperature range
	Do not resterilize		Keep dry
<b>REF</b>	Catalogue or model number	<b>LOT</b>	Batch code
	Do not use if package is damaged		Date of manufacture
<b>MD</b>	Medical Device	<b>UDI</b>	Unique Device Identification
<b>Rx ONLY</b>	The medical device according to 21 CFR 801.109(b)(1) may only be prescribed by a physician in accordance with the law (USA)		

© The NGMedical brands and registered trademarks are protected by national and international rights in Germany, the European Union and other regions and countries around the world.

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NGMedial GmbH