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PLIF Implantate

DE

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EN

■ INTENDED USE

BEE® PLIF cages and NGM WAVE® PLIF cages are implants for intervertebral, posterior, lumbar fusion in the segments L2-S1 among patients with a mature skeletal system.

■ PRODUCT DESCRIPTION

BEE® PLIF cages and NGM WAVE® PLIF cages have been specially adapted to suit the local anatomy to ensure the best possible surgical outcome. BEE® PLIF cages are made of titanium alloy in accordance with ASTM F3001 and are produced in an additive manufacturing process. NGM WAVE® PLIF cages are made of PEEK-OPTIMA™ Natural and tantalum (pins) in accordance with ASTM F560. The implants are available in various heights, lengths and lordosis angles. The cranial and caudal side of the implant is rough or fitted with teeth to achieve a good anchoring of the implant.

■ SCOPE

These instructions for use are valid for all BEE® PLIF and NGM WAVE® PLIF implants.

■ APPROVED INSTRUMENTS

Only the instruments described in the BEE® PLIF and NGM WAVE® PLIF surgical techniques may be used in direct connection with the implant. The surgical technique is available from NGMedical on request.

■ INDICATIONS

BEE® PLIF and NGM WAVE® PLIF are intended for the following indications:

- degenerative disorders of the spine including degenerative disc disease, vertebral canal and foraminal stenosis,
- lumbar spine instability in one or several segments,
- pseudoarthrosis,
- spondylolisthesis, including degenerative spondylolisthesis,
- spondylodiscitis.

■ CONTRAINDICATIONS

BEE® PLIF and NGM WAVE® PLIF must not be implanted in the following patients or health conditions:

- obesity from grade II (BMI value ≥ 35),
- active systemic infection,
- allergy to the implant material (see the imprint on the packaging).
- endocrine or metabolic illness in the patient's history which influences the bone and mineral metabolism,
- health conditions which are opposed to the potential benefit of a spinal operation,
- the patient has received medication influencing bone and mineral metabolism less than two weeks before the planned date of operation,
- osteoporosis,
- the patient is unwilling or unable to follow the post-operative instructions,
- pregnancy,
- vertebral fracture(s),
- spinal metastases (exception: additional surgical stabilization measures).

■ COMPLICATIONS AND POSSIBLE ADVERSE EFFECTS

Before the operation, the possible adverse effects and possibly required further surgery should be explained to the patient:

- leakage of cerebrospinal fluid,
- spinal curvature with kyphotic or scoliotic deformity, loss of correction,
- (accelerated) adjacent segment degeneration,
- breakage/failure of system components (loosening, degradation or breakage of components),
- permanent disability,
- dislocation/migration of the implant,
- subsidence of the implant into the vertebral body end plates causing reduced height of the intervertebral space,
- failed fusion with pseudoarthrosis,
- foreign body reactions to components, their debris or their fragments,
- infection,
- loss of bone density due to stress shielding,
- possibly required further surgery,
- nerve damage and/or spinal cord damage including the loss of neurological functions (sensory and/or motor deficiencies such as hypoaesthesia, anesthesia, dysesthesia, hyperaesthesia, paresthesia, partial or complete muscular paralysis, lack of reflexes, increased or reduced muscular tonus, spasticity),
- osteoporotic vertebral fracture,
- spinal cord injuries with partial or complete paraplegia,
- death,
- visceral or vascular injuries,
- spinal infections (e.g. spondylitis, discitis, myelo-meningitis)
- wound healing disorders.

■ WARNINGS

The safety and effectiveness of the BEE® PLIF and NGM WAVE® PLIF implants have been proven only for the spinal disorders stated in the indications section. As with all operations, the use of this system involves the potential risk of death. Other potential risks which may require further surgery are as follows:

- allergic reaction or sensitivity to the implant material (see the imprint on the packaging).
- breakage/failure of system components,
- vascular injuries,
- pseudoarthrosis,
- damage to the spinal cord and/or nerves,
- loss of anchorage.

By implantation of metals and alloys into the human body, these are exposed to aggressive chemical milieus of salts, acids and alkalis which might lead to corrosion. If different metals are positioned in close vicinity, galvanic effects may accelerate the corrosion process. In particular, BEE® PLIF and NGM WAVE® PLIF must not be used with components (exception: bone substitute material or fillers) or instruments from other manufacturers. In general, implanted material may cause immune reactions or foreign body reactions as well as chronic inflammatory changes.

When removing the implant from the packaging, it must be ensured that the nominal height of the implant matches the size stated on the packaging labeling.

■ PRECAUTIONS

The surgeon must be completely familiar with BEE® PLIF and NGM WAVE® PLIF, the required instruments and the general surgical technique. The surgeon must use the instruments which are specifically approved by NGMedical for use with BEE® PLIF and NGM WAVE®. The surgeon must always ensure that the adequate implant size is used and that the implant is fitted into the anatomically correct position. Overdistraction must be avoided. It must be observed that the implant is aligned correctly in the cranial-caudal as well as the ventral-dorsal plane. The steps of the operation can also be read in the implant-specific surgical technique.

Documentation of the used implants is necessary for each patient with REF, LOT and UDI so that the level of traceability required by law is guaranteed. The implants are intended for single use. An implant must not be reused once explanted.

The BEE® PLIF and NGM WAVE® PLIF implants have not been tested with regard to their safety and compatibility with MRI. They have also not been tested for heating or migration in an MRI environment.

■ PACKAGING

BEE® PLIF and NGM WAVE® PLIF implants are generally supplied in double sterile packaging. Only use the implants if the packaging is intact. The storage conditions must guarantee the intact nature of the packaging and therefore also the implants. Defective implants must not be used. Each packaging has a documentation label for you to use and an implantation ID. The packaging and labels contain all the required information. There, you can find the product number REF, LOT, description, method of sterilization, quantity, expiration date and manufacturer, as well as scan the UDI code. All BEE® PLIF and NGM WAVE® PLIF implants are sterilized by way of gamma irradiation sterilization.

■ COMPLAINTS


Every specialist (e.g. doctor who uses the product) who has a complaint or who is dissatisfied with the quality, labeling, reliability, safety, effectiveness and/or performance of the product should inform NGMedical or, where applicable, the distributor. In the case of a severe incident or the risk of a severe incident which may lead or has led to a severe health impairment or the death of the patient, a user or third-party, NGMedical (or the retailer) must be informed immediately either verbally or in writing. All complaints must be made, stating the product name, REF, LOT, UDI and description as well as the date of the incident. The reporting party must state their name and contact information, and describe the incident in as much detail as possible. The legal reporting obligations must also be observed!

■ FURTHER INFORMATION

For further information, please contact:

NGMedical GmbH
Talsperrenblick 5
66620 Nonnweiler
GERMANY
mail@ngmedical.de
www.ngmedical.de
+49 (0) 6873 99997-0

	Consult instructions for use	STERILE R	Sterilized using irradiation
	Use by		Manufacturer
	Do not re-use		Storage temperature range
	Do not re-sterilize		Keep dry
REF	Catalogue or model number	LOT	Batch code
	Do not use if package is damaged		Date of manufacture
MD	Medical Device	UDI	Unique Device Identification
Rx ONLY	The medical device according to CFR 820.3 (I) may only be prescribed by a physician in accordance with the law (USA)		

© The NGMedical brands and registered trademarks are protected by national and international rights in Germany, the European Union and other regions and countries around the world.

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