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BEE® HA Implantate

REF X0005 LOT 2025-09-24

DE

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EN-NOT FOR US

■ INTENDED PURPOSE

BEE® HA Cages were developed as an intercorporeal implant for the anterior cervical spondylosis (C3-C7) among patients with a mature skeletal system.

■ PRODUCT DESCRIPTION

BEE® HA Cages were especially adapted to the local anatomy to secure the surgical outcome at the best. The caudal side is flat, the cranial side is domed and the implant is formed conically from anterior to posterior. In the lateral view, the implant has a slightly lordotic form and is available in 4°, 8° and 12° lordosis. The cage disposes of 4 pins on both bone facing surfaces, which serve the primary fixation. BEE® HA implants are made of PEEK-OPTIMA® HA ENHANCED and pins made out of titanium alloy (Ti6Al4V_ELI) according to ASTM F136.

It is not allowed to use BEE® HA implants in contact with components of other manufacturers with the exception of bone graft substitutes or fillers.

■ SCOPE

These Instructions for Use are valid for all cervical BEE® HA implants made of PEEK-OPTIMA® HA ENHANCED with pins made of titanium alloy (Ti6Al4V_ELI).

■ INTENDED USERS

These instructions for use alone do not provide sufficient information for immediate use of the system. Training in the use of the products by an experienced and trained surgeon and/or by NGMedical (or the distributor) is strongly recommended.

The procedure must be performed in accordance with the instructions for use and the surgical technique. The surgeon is responsible for ensuring that the operation is performed correctly. It is strongly recommended that the procedure is only performed by surgeons who are qualified and experienced in spinal surgery, are aware of the general risks associated with spinal surgery and are familiar with the product-specific surgical techniques.

■ APPROVED INSTRUMENTS

Only the instruments described in the cervical surgical technique (REF X1003 (German), X1004 (English) or X1008 (US)), may be used in direct connection with the implant.

■ EXPECTED CLINICAL BENEFIT

The cervical cage BEE® HA is intended for use in accordance with its intended purpose and instruction for use for cervical fusion of adjacent vertebral bodies. The expected clinical benefits of surgery with BEE® HA, when used as intended, include a low complication rate (especially subsidence rate) and good osseointegration (fusion) due to optimized design features such as shape and anchoring, resulting in satisfactory radiological results. Furthermore, improvement in pain and quality of life for the patient is expected. The clinical benefits are supported by preclinical data and study data.

The Summary of Safety and Clinical Performance (SSCP) can be found under the following link: www.ngmedical.de/sscp/bee_ha

Access to the public website of the Eudamed database: <http://ec.europa.eu/tools/eudamed>.

■ INDICATIONS

BEE® HA is intended to be used for the following indications:

- Discopathy,
- Disc hernia with radiculopathy and / or myelopathy,
- Stenosis of the foramen and spinal canal with or without cervical radiculopathy/ myelopathy.

■ CONTRAINDICATIONS

It is not allowed to implant BEE® HA in patients with the following health status:

- Bone mineral density with T-score \leq -2.5 as determined by spine DXA if male \geq 60 years of age or female \geq 50 years of age (Exception: additional operational stabilization measures),
- Active systemic infection or infection at the surgical site,
- Sustained osteoporotic fracture/ osteopenia of the spine, hip or wrist (Exception: additional operational stabilization measures),
- Medications (e.g., methotrexate, alendronate) that interfere with bone and mineral metabolism within 2 weeks of the planned date of the index surgery,
- Any medical or surgical condition precluding the potential benefit of spinal surgery,
- History of endocrine or metabolic disorders (e.g., Paget's disease) known to affect bone and mineral metabolism,
- Spinal metastases (Exception: additional operational stabilization measures),
- Known allergy to polyetheretherketone (PEEK), hydroxyapatite (HA) or titanium alloy (Ti6Al4V_ELI),
- Pregnancy,
- Axial neck pain as the solitary symptom,
- Advanced cervical anatomical deformity (e.g., ankylosing spondylitis, scoliosis) at the operative site,
- The patient is not willing or able to follow the postoperative instructions.

■ COMPLICATIONS AND POSSIBLE UNDESIRE OCCURENCES

Prior to the surgery, patients should be explained the possible undesired occurrences and the possible further surgeries which might be necessary for the correction of these:

- Breakage/ failure of system components (loosening, collapse or breaking of components),
- Dislocation/ migration of the implant,
- Possible wound healing disorder and/or wound complications,
- Infection,
- Neural damages and/or spinal cord damages including loss of neurological functions (sensory and/or motor deficits such as hypoaesthesia, anaesthesia, dysaesthesia, hyperaesthesia, paraesthesia, muscle paresis or paralysis, reflex deficiency, increased or decreased muscle tone, spasticity),
- Spinal cord injuries with partial or complete transectional lesion,
- Leakage of the cerebrospinal fluid,
- Spinal fluid infection (myelo-meningitis),
- Trachea and/or oesophagus injury,
- Vocal cord paresis,
- Hoarseness,
- Foreign body reaction to components or fragments/debris,
- Visceral or vascular injuries (Carotid arteries, jugular veins),
- Modification of the spinal curvature with kyphotic or scoliotic deformity, loss of correction,
- Loss of bone density by stress-shielding,
- Need of further surgeries,
- Failed fusion with pseudarthrosis,
- Subsidence of the implant into the endplates with loss of intervertebral space height,
- Accelerated degeneration of adjacent segments (Adjacent Segment Disease),
- Permanent disability,
- Osteolysis, bone loss or fracture, bone resorption above or below the level of the surgical segment,
- Dysphagia,
- Death.

■ WARNINGS

The safety and efficacy of the BEE® HA implants has only been proven for the spinal diseases listed in the paragraph „indications“. There is the potential risk of death when using this system. Further potential risks, which might lead to further surgeries, are:

- Breakage/ Failure of system components,
- Loss of Anchorage,
- Pseudarthrosis,
- Allergic reaction or sensitivity on polyetheretherketone (PEEK), hydroxyapatite (HA) or titanium alloy (Ti6Al4V_ELI),
- Nerve damage and/or spinal cord damage,
- Vascular lesion.

By implantation of metals and alloys into the human body, these are exposed to aggressive chemical milieu of electrolytes, acids and alkalis, which might lead to corrosion. Particularly, it is not allowed to use the BEE® HA Cage with components (except for bone graft substitutes or fillers) or surgical techniques of other manufacturers. Generally, implanted materials may cause immune responses or foreign-body reactions with chronic inflammatory changes. It must be ensured that the size of the cage itself matches with the size indicated on the packaging label, when taking the device of its packaging box.

■ PRECAUTIONS

The surgeon must be perfectly familiar with BEE® HA, the required instruments as well as the general surgical technique. The surgeon must exclusively use the instruments especially approved by NGMedical for BEE® HA. The surgeon always has to ensure that the optimum implant size is used and that there is an anatomically correct implant placement. Trial implants enable a simple and safe size selection. Over-distraction is to be avoided. The surgeon has to ensure the correct implant placement in the cranial, caudal, lateral, anterior and posterior position. The general surgical procedure is described in detail in the surgical technique provided by NGMedical (REF X1003 (German), X1004 (English) or X1008 (US)). Implanted devices are to be documented patient-related with REF, LOT and UDI-

Code, so that the legally required traceability is guaranteed. The implants are intended for single use. **An explanted implant must never be reused.**

■ MRI COMPATIBILITY

The effects of magnetic resonance (imaging) on BEE® HA implants have not been tested. Based on the widespread use of titanium alloy implants (Ti6Al4V_ELI), BEE® HA implants can be categorized as having limited MRI suitability. Patients with these implants can be scanned under 1.5 Tesla and 3.0 Tesla static magnetic field strength. A maximum specific whole-body absorption rate (SAR) of 1.75 W/kg at 1.5 Tesla and 3.0 Tesla for a scan duration of 15 minutes should not be exceeded.

■ PACKAGING

BEE® HA implants are generally supplied in double sterile packaging. Please only use implants if the packaging box is undamaged. Storage conditions must guarantee the integrity of the packaging box and the implants respectively. Defective implants or other components are not to be used. Each packaging box contains documentation labels for free use and an implant card. The packaging box and the labels contain all information required. They indicate the REF, LOT, description, sterilization method, quantity, expiration date, manufacturing date, manufacturer and UDI-Code. All BEE® HA implants are sterilized using gamma radiation or ethylene oxide gas. The type of sterilization is indicated on each packaging box.

■ COMPLAINTS

Any healthcare professional (e. g., a surgeon using the product), who has a complaint or who is dissatisfied with the quality, identification, reliability, safety, efficacy and/or performance of the BEE® HA products, should notify NGMedical or, where applicable, the distributor. In the case of a severe incident or risk of a severe incident, liable to result in or to have resulted in the death or severe deterioration in the state of health of a patient or user or other persons, NGMedical (or the distributor) and the competent authority must be informed as soon as possible in writing or verbal. All complaints must include the product name, REF, LOT, UDI and the description as well as all data of the incident. The reporting person should state his name, contact data and a detailed description of the incident. The legal reporting obligations are also to be observed!

■ DISPOSAL

For disposal of NGMedical implants, please comply with the specific hospital regulations or return the processed, sterilized and labeled implants to NGMedical.

■ EXPECTED LIFETIME


The expected lifetime of BEE® HA is 2.5 years after implantation under non-fused conditions. After successful fusion, the expected lifetime corresponds to the remaining lifetime of the patient. If the fusion is not successful, a revision may become necessary.

■ FURTHER INFORMATION

For further information please contact

NGMedical GmbH
Trierer Straße 65
66620 Nonnweiler
GERMANY
mail@ngmedical.de
www.ngmedical.de
+49 (0) 6873 99997-0

Symbol	Standard	Ref no	Title	Meaning
		5.4.3	Consult instructions for use	Indicates the need for the user to consult the instructions for use.
		5.1.4	Use by	Indicates the date after which the medical device is not to be used.
		5.4.2	Do not re-use	Indicates a medical device that is intended for one use, or for use on a single patient during a single procedure.
		5.2.6	Do not re-sterilize	Indicates a medical device that is not to be re-sterilized.
		5.1.6	Catalogue or model number	Indicates the manufacturer's catalogue number so that the medical device can be identified.
		5.2.8	Do not use if package is damaged	Indicates a medical device that should not be used if the package has been damaged or opened.
	ISO 15223-1- Symbols to be used with information to be supplied by the manufacturer - Part 1: General requirements	5.2.4	Sterilized using irradiation	Indicates a medical device that has been sterilized using irradiation.
		5.2.3	Sterilized with ethylene oxide	Indicates a medical device that was sterilized with ethylene oxide
		5.1.1	Manufacturer	Indicates the medical device manufacturer.
		5.3.7	Storage temperature range	Indicates the temperature limits to which the medical device can be safely exposed.
		5.3.4	Keep dry	Indicates a medical device that needs to be protected from moisture.
		5.1.5	Batch code	Indicates the manufacturer's batch code so that the batch or lot can be identified.
		5.1.3	Date of manufacture	Indicates the date of manufacture the medical device
		5.7.7	Medical Device	Indicates the item is a medical device
		5.7.10	Unique Device Identification	Indicates a carrier that contains Unique Device Identifier information
			Double sterile barrier system	
			European conformity mark with Notified Body identification number	
			The medical device according to 21 CFR 801.109(b)(1) may only be prescribed by a physician in accordance with the law (USA)	

INSTRUCTIONS FOR USE

BEE® HA Implants

REF

X0005

LOT

2025-09-24

	Conditionally MRI compatible
---	------------------------------

© The NGMedical brands and registered trademarks are protected by national and international rights in Germany, the European Union and other regions and countries worldwide.

INSTRUCTIONS FOR USE

BEE® HA Implants

REF X0005 LOT 2025-09-24

EN | FOR US ONLY

■ CAUTION

Federal law (USA) restricts these devices to sale by or on the order of a physician. These devices should be implanted only by a physician who is fully trained with the devices, intended use, instrumentation and with knowledge of the surgical techniques required. Contact your representative for surgical technique.

■ INTENDED PURPOSE

BEE® HA Cages were developed as an intercorporeal implant for the anterior cervical spondylosis among patients with a mature skeletal system.

■ PRODUCT DESCRIPTION

BEE® HA Cages were especially adapted to the local anatomy to secure the surgical outcome at the best. The caudal side is flat, the cranial side is domed and the implant is formed conically from anterior to posterior. In the lateral view, the implant has a slightly lordotic form and is available in 4°, 8° and 12° lordosis. The cage disposes of 4 pins on both bone facing surfaces. BEE® HA implants are made of PEEK-OPTIMA® HA ENHANCED and pins made out of titanium alloy (Ti6Al4V_ELI) according to ASTM F136.

It is not allowed to use BEE® HA implants in contact with components of other manufacturers with the exception of autogenous and/or allogeneic bone graft comprised of cancellous and/or corticocancellous bone graft.

■ SCOPE

These Instructions for Use are valid for all cervical BEE® HA implants made of PEEK-OPTIMA® HA ENHANCED with pins made of titanium alloy (Ti6Al4V_ELI) according to ASTM F136.

■ APPROVED INSTRUMENTS

Only the instruments described in the BEE® HA surgical technique (X1008 (US)), may be used in direct connection with the implant. The surgical technique is available on request from NGMedical.

■ INDICATIONS

BEE® HA Cages are intended for intervertebral body fusion devices in skeletally mature patients for the treatment of cervical disc degeneration and/or cervical spinal instability as confirmed by imaging studies (radiographs, CT, MRI) that result in radiculopathy, myelopathy and/or pain at one or more contiguous levels from C2-T1. These patients should have had at least six weeks of nonoperative treatment. BEE® HA Cages are to be used with autogenous and/or allogeneic bone graft comprised of cancellous and/or corticocancellous bone graft to facilitate fusion and in combination with supplemental fixation indicated for cervical fusion procedures.

■ CONTRAINDICATIONS

It is not allowed to implant BEE® HA in patients with the following health status:

- Bone mineral density with T-score ≤ -2.5 as determined by spine DXA if male ≥ 60 years of age or female ≥ 50 years of age (Exception: additional operational stabilization measures),
- Active systemic infection or infection at the site,
- Sustained osteoporotic fracture/ osteopenia of the spine, hip or wrist (Exception: additional operational stabilization measures),
- Medications (e.g., methotrexate, alendronate) that interfere with bone and mineral metabolism within 2 weeks of the planned date of the index surgery,
- Any medical or surgical condition precluding the potential benefit of spinal surgery,
- History of endocrine or metabolic disorders (e.g., Paget's disease) known to affect bone and mineral metabolism,
- Spinal metastases (Exception: additional operational stabilization measures),
- Known allergy to polyetheretherketone (PEEK), hydroxyapatite (HA) or titanium alloy (Ti6Al4V_ELI),
- Pregnancy,
- Axial neck pain as the solitary symptom,
- Advanced cervical anatomical deformity (e.g., ankylosing spondylitis, scoliosis) at the operative site,
- The patient is not willing or able to follow the postoperative instructions,
- Active infections or high risk of infection,
- Signs of local inflammation,
- Fever or leukocytosis,
- Morbid obesity,
- Mental illness,
- Any patient having inadequate tissue coverage over the operative site,
- Use of this type of surgical implant surgery in children or pediatric patients presents particular risks because of bone growth or physical movement. Subsequent re-intervention may be required,
- Infants with a known hereditary or acquired bone friability or calcification problem, or those with a very short life expectancy, should not be considered for this type of surgery,
- Substance abuse or senility that precludes the patient from following post operative precautions to prevent implant failure,
- Any patient in which implant utilization would interfere with anatomical structures or expected physiological performance.

■ COMPLICATIONS AND POSSIBLE UNDESIRE OCCURENCES

Prior to the surgery, patients should be explained the possible undesired occurrences and the possible further surgeries which might be necessary for the correction of these:

- Breakage/ failure of system components (loosening, collapse or breaking of components),
- Dislocation/ migration of the implant is higher without the use of supplemental fixation,
- Possible degeneration of the skin and/or wound complications,
- Infection,
- Nerval damages and/or spinal cord damages including loss of neurological functions (sense or mobility paralysis), dysaesthesia, hyperaesthesia, paraesthesia, radiculopathy, reflex deficiency,
- Spinal cord injuries, outflow of the cerebrospinal fluid,
- Trachea and/or esophagus injury,
- Foreign body reaction (allergically based) to components or fragments/debris,
- Visceral or vascular injuries,
- Modification of the spinal curvature. Loss of correction, height and/or repositioning,
- Loss of bone density by stress-shielding,
- Need of further surgeries,
- Failed fusion,
- Subsidence of the implant into the endplates,
- Adjacent Segment Disease,
- Osteolysis, bone loss or fracture, bone resorption above or below the level of the surgical segment,
- Dysphagia,
- Death.

■ WARNINGS

The safety and efficacy of the BEE® HA implants have only been proven for the spinal diseases listed in the paragraph „indications“. There is the potential risk of death when using this system. Further potential risks, which might lead to further surgery, are:

- Breakage/Failure of system components,
- Loss of Anchorage,
- Pseudarthrosis,
- Allergic reaction or sensitivity to polyetheretherketone (PEEK), hydroxyapatite (HA) or titanium alloy (Ti6Al4V_ELI),
- Nerve damage and/or spinal cord damage,
- Vascular lesion.

By the implantation of metals and alloys into the human body, these are exposed to aggressive chemical milieu of salts, acids and alkalis, which might lead to corrosion. Particularly, it is not allowed to use the BEE® HA Cage with components or surgical techniques of other manufacturers. It must be ensured that the size of the cage itself matches with the size indicated on the packaging label, when taking the cage out of its packaging box.

■ PRECAUTIONS

The surgeon must be perfectly familiar with BEE® HA, the required instruments as well as the general surgical technique. The surgeon must exclusively use the instruments especially approved by NGMedical for BEE® HA. The surgeon always has to ensure that the optimum implant size is used and that there is an anatomically correct implant placement. Trial implants enable a simple and safe size selection. Over distraction is to be avoided. The surgeon must ensure that there is a correct implant placement in the cranial and caudal as well as in the anterior and posterior position. The general surgical procedure is described in detail in the surgical technique (REF X1008). This document can be obtained from NGMedical or the representative.

Single use only. Implants are provided as single use only implants and are not to be reused or implanted regardless of an apparent undamaged condition.

Please document the used implants patient-related with REF and LOT and UDI-Code, so that the legally required traceability is guaranteed.

Based on the fatigue testing results, the physician/surgeon should consider the levels of implantation, patient weight, patient activity level, other patient conditions, etc. which may impact on the performance of the system.

■ MRI COMPATIBILITY

BEE® HA has not been evaluated for safety and compatibility in the MR environment. It has not been tested for heating migration, or image artifact in the MR environment. The safety of the BEE® HA system in the MR environment is unknown. Scanning a patient who has this device may result in patient injury.

■ PACKAGING

BEE® HA implants are generally supplied in double sterile packaging. Please only use implants if the packaging box is undamaged. Storage conditions must guarantee the integrity of the packaging box and the implants respectively. Defective implants or other components are not to be used. Each packaging box contains documentation labels for free use and an implant card. The packaging box and the labels contain all information required. They indicate the REF, LOT, description, sterilization method, quantity, expiration date, manufacturing date, manufacturer and UDI-Code. All BEE HA implants are sterilized by gamma radiation.

■ COMPLAINTS

Any healthcare professional (e.g., a surgeon using the product), who has a complaint or who is dissatisfied with the quality, identification, reliability, safety, efficacy and/or performance of the BEE® HA products, should notify NGMedical or, where applicable, the distributor. In the case of a severe incident or risk of a severe incident, liable to result in or to have resulted in the death or severe deterioration in the state of health of a patient or user or other persons, NGMedical (or the distributor) has to be informed as soon as possible in writing or verbal. All complaints must include the product name, REF, LOT, UDI and the description as well as all data of the incident. The reporting person should state his name, contact data and a detailed description of the incident. The legal reporting obligations are also to be observed!

■ DISPOSAL

For disposal of NGMedical implants, please comply with the specific hospital regulations or return the processed, sterilized and labeled implants to NGMedical.

■ EXPECTED LIFETIME


The expected lifetime of BEE® HA is 2.5 years after implantation under non-fused conditions. After successful fusion, the expected lifetime corresponds to the remaining lifetime of the patient. If the fusion is not successful, a revision may become necessary.

■ FURTHER INFORMATION

For further information please contact

NGMedical GmbH
Trierer Straße 65
66620 Nonnweiler
GERMANY
mail@ngmedical.de
www.ngmedical.de
+49 (0) 6873 99997-0

The Summary of Safety and Clinical Performance (SSCP) can be found under the following link: www.ngmedical.de/sscp/bee_ha.

Access to the public website of the Eudamed database: <http://ec.europa.eu/tools/eudamed>.

Symbol	Standard	Ref no	Title	Meaning
		5.4.3	Consult instructions for use	Indicates the need for the user to consult the instructions for use.
		5.1.4	Use by	Indicates the date after which the medical device is not to be used.
		5.4.2	Do not re-use	Indicates a medical device that is intended for one use, or for use on a single patient during a single procedure.
		5.2.6	Do not re-sterilize	Indicates a medical device that is not to be re-sterilized.
		5.1.6	Catalogue or model number	Indicates the manufacturer's catalogue number so that the medical device can be identified.
		5.2.8	Do not use if package is damaged	Indicates a medical device that should not be used if the package has been damaged or opened.
	ISO 15223-1- Symbols to be used with information to be supplied by the manufacturer - Part 1: General requirements	5.2.4	Sterilized using irradiation	Indicates a medical device that has been sterilized using irradiation.
		5.1.1	Manufacturer	Indicates the medical device manufacturer.
		5.3.7	Storage temperature range	Indicates the temperature limits to which the medical device can be safely exposed.
		5.3.4	Keep dry	Indicates a medical device that needs to be protected from moisture.
		5.1.5	Batch code	Indicates the manufacturer's batch code so that the batch or lot can be identified.
		5.1.3	Date of manufacture	Indicates the date of manufacture the medical device
		5.7.7	Medical Device	Indicates the item is a medical device
		5.7.10	Unique Device Identification	Indicates a carrier that contains Unique Device Identifier information
		Double sterile barrier system		
		European conformity mark with Notified Body identification number		
		The medical device according to 21 CFR 801.109(b)(1) may only be prescribed by a physician in accordance with the law (USA)		
		Conditionally MRI compatible		

© The NGMedical brands and registered trademarks are protected by national and international rights in Germany, the European Union and other regions and countries worldwide.

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