Aufbereitungsanweisung Instrumente / Reprocessing Instructions Instruments / Instructions de préparation d'instruments / Istruzioni di purificazione degli strumenti / Instruções de reprocessamento dos instrumentos / Instrucciones de tratamiento de los instrumentos / 医疗器械准备说明 / तैयारी के निर्देश उपकरण / 器材の再処理説明書 / Инструкция по обработке инструментов / Οδηγίες προετοιμασίας εργαλείων

REF X0200 LOT 2024-04-12

DE

1. ALLGEMEINE GRUNDLAGEN

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EN

1. GENERAL BASICS

In the following, a validated procedure is presented, which describes the cleaning, disinfection, and sterilization of the reprocessable instruments (critical B, according to Spaulding and RKL) from NGMedical. Prior to each use, the process must be performed completely. In case of any deviation from the described process, NGMedical does not take over responsibility for the sterility and cleanliness of the instruments.

Additionally, please respect the legal regulations valid in your country as well as the hygiene rules of the hospital. This especially applies to the various specifications regarding an effective inactivation of prions.

Prior to each use, all instruments must be cleaned, disinfected, and sterilized. This is especially applicable for the initial use after delivery, as all instruments are delivered unsterile, but are free from production residues (cleanliness and disinfection after removal of the protective transport packaging; sterilization after positioning in sterilization tray). A proper cleaning and disinfection are an essential prerequisite for an effective sterilization of the instruments.

Within your responsibility for the fulfilment of all steps of this instruction, please see to it that only sufficiently device- and product-specific validated procedures are used for the cleaning/disinfection and sterilization, that the used equipment (ultrasound, cleaning and disinfection device, sterilization device) are maintained and inspected regularly and that the validated parameters are observed for each cycle.

When using the instruments, please already pay attention to collect contaminated instruments separately and not placed in the instrument tray again to avoid a higher contamination of the equipped instrument tray.

It should also be noted that there may also be smaller instrument trays in the instrument tray. Please also remove the unused instruments from these trays for cleaning/disinfection.

The time between using the instruments and cleaning must not exceed two (2) hours or the soiling must not be allowed to dry. Clean/disinfect the soiled instruments, then sort them back into the instrument tray and sterilise the fully loaded instrument tray.

Some instruments require special aspects. Please see chapter "Special Instructions".

2. CLEANING AND DISINFECTION

Initial treatment at the site of use

At the site of use, please proceed as follows:

- Remove pulp residues by wiping or rinsing with distilled water immediately after use
- Put the roughly cleaned instruments in an appropriate storage tray (wet storage) and, if necessary, mark damaged or defective instruments
- Use a cleaning and/or disinfectant solution. The cleaning and/or disinfectant solution must be aldehyde-free (otherwise blood rests might not be removed), have a proven effectiveness (e. g. DGHM – or FDA-approval or CE-approval respectively), must be applicable for the disinfection of the instruments and be compatible with the instruments (see chapter "Material Resistance").
- For the manual removal of residues, only use a soft brush or a clean, soft tissue, which are only used for this purpose, however, never use metal brushes or steel wool.
- Movable instruments need to be moved several times during pre-cleaning (e. g. joints)
- If applicable: Rinse all lumina of the instruments by using a single use syringe (min. volume 10 ml)

Please note that the disinfectant solution, which is used during pre-treatment, only serves as a personal protection and cannot replace the following disinfection step which is carried out after pre-cleaning. All instruments, including the instruments, which have not been used, must be reprocessed. Also, it is not allowed to leave out any step of reprocessing in case of not used instruments. Additionally, a pre-cleaning has to be affected in an ultrasonic bath.

Machine Cleaning and Disinfection Device

For ultrasonic cleaning, please take care

- That the instruments do not touch each other
- That the instruments ideally are placed into the ultrasonic bath preventing rinse residues
- That all instruments are - interiorly and exteriorly – completely covered with liquid.
- That the detergent agent must be freshly prepared every working day as a minimum and must be changed immediately in case of visible contamination
- That the detergent agent is prepared according to the instructions of the manufacturer
- That the purification bathes are cleaned every working day thoroughly mechanically and disinfected in order to avoid biofilm formation

When choosing the cleaning and disinfection device please take care

- That the cleaning and disinfection device principally disposes of a proven effectiveness (e.g. DGHM – or FDA-approval, or CE-approval respectively, according to DIN EN ISO 15883).
- That a proven program for thermal disinfection (min. 5 min at 90°C or A0-value >3000) is applied (in case of chemical disinfection risk of disinfection solution residues on the instruments)
- That the used program is applicable for all instruments and includes a sufficient number of rinsing cycles
- That for rinsing, ideally sterile or low-germ (max. 10 germs/ml) as well as low endotoxin (max 0.25 endotoxin units/ml) water (e. g. highly purified water) is used

² longer sterilization period for prion inactivation (e. g. 18 min at 134°C (273,5°F) in Germany). Please follow your national rules.

- That the drying air is filtered and
- That the cleaning and disinfection device is maintained and inspected regularly

For the selection of the detergent systems please take care

- that it is principally applicable for cleaning instruments out of metal and plastics
- that the chemicals used are compatible with the instruments (see chapter "material resistance") and
- that a cleaning agent with a pH>10 to 11 is used (e. g. neodisher MediClean forte of Dr. Weigert)

It is essential that the concentrations of the detergent and disinfection agents, indicated by the manufacturer, are strictly observed.

Process

1. Disassemble the instruments (please observe chapter "Special Instructions").
2. Place the disassembled, in case of instruments with joints, the opened instruments, into the ultrasonic bath. Please see to it that the instruments do not touch each other and that there are no acoustical shadows.
3. All instruments must be completely sprinkled interiorly and exteriorly, cavities must be filled free of air bubbles (please insert slightly angled)
4. Start the ultrasonic cleaning of the instruments with the following parameter:
 - a. Frequency 32 to 50 kHz
 - b. Temperature max. 40°C (104°F)
 - c. 10 minutes
5. During pre-cleaning in the ultrasonic bath, it is not allowed to insert further instruments or subjects
6. Take out the instruments when the program of the ultrasonic bath is finished
7. Inspect the instruments regarding visible modifications and contamination
8. Put all the disassembled instruments into the cleaning and disinfection device. Please see to it that the instruments do not touch each other.
9. Adapt all lumina of the instruments to the rinse connection of the cleaning and disinfection device by using an appropriate rinse adapter.
10. Start the program with the following parameter:
 - a) 1 min. pre-cleaning with cold (<40°C (104°F)) water (drinking water quality according to EG guideline 98/83/EG)
 - b) 10 min. cleaning with e. g. Mediclean 0.5% (Dr. Weigert) at 55°C (131°F)
 - c) 2 min neutralization with neodisher Z 0.1% (Dr. Weigert)
 - d) 2 min rinsing with cold (<40°C (104°F)) completely desalinated water (demineralized water)
 - e) Thermal disinfection, min. 5 min at 90°C (194°F) or A0-values-3000
11. After program finish of the cleaning and disinfection device, take out the instruments
12. Inspect and pack the instruments as soon as possible after removal (see chapter "inspection", "maintenance" and "packaging", eventually after additional drying, at a clean place).

Some instruments require special aspects. Please see chapter "Special Instructions".

The fundamental suitability of an effective automated cleaning and disinfection was demonstrated by an independent accredited test laboratory by application of the RDG Miele PG 8535 and the specified procedure described above.

3. INSPECTION

Inspect all instruments after pre-cleaning and cleaning/disinfection regarding the following characteristics:

- Corrosion
- Damaged surfaces
- Split-offs
- Contamination
- Fissures and sharp edges
- Deforming

Sort out the instruments indicating one or several of the a.m. characteristics. Instruments which are still contaminated must be cleaned and disinfected again. If the instrument is still contaminated after a repeated cleaning/disinfection, it must be removed from circulation.

4. MAINTENANCE

For instruments with complex functions (e.g. ratchets, couplings), it is necessary to lubricate movable components and running surfaces after each reprocessing cycle. Each insufficient maintenance, e.g. by not lubricating or not sufficiently lubricating, might lead to a damage of the instrument. When using instrument oils, please see to it that only instrument oils (white oils) are applied, which are - under consideration of the maximum applied sterilization temperature - approved for steam sterilization and dispose of a proven biocompatibility (e.g. Steril from Aesculap). Put instrument oil punctually on movable parts (joints and sliding surfaces) and spread the oil by moving the parts.

Some instruments require special aspects. Please see chapter "Special Instructions".

5. PACKAGING

Sort the cleaned and disinfected instruments into the corresponding sterilization tray (see chapter "special indications"). The instruments are not allowed to touch each other in the tray.

Please pack the instruments and the sterilizations trays in sterilization containers, which meet the following requirements:

- According to EN ISO 11607/ ANSI/AAMI TIR 11607
- Applicable for steam sterilization (temperature resistance up to min. 141°C, (286°F) sufficient steam permeability)
- Sufficient protection of the instruments and sterilization packaging from damages
- Maintained regularly according to the manufacturer's indications (sterilization container)

The sets must be packed in accordingly sealed sterilization trays/containers. Sterilization cloths do not offer a sufficient protection.

The fundamental suitability of an effective steam sterilization was demonstrated by an independent accredited test laboratory by application of a sterilization container Ermis MedTech GmbH. The procedure described here, taking into account the worst case parameters, was taken into account.

6. STERILIZATION

Only the following mentioned sterilization methods are to be used for the sterilization; other sterilization method are not permitted:

- Steam sterilization
- Fractionated vacuum method (minimum three (3) cycles with sufficient product drying)
- Steam sterilizer according to EN 13060/ ANSI/AAMI ST79 or EN 285/ ANSI/AAMI ST79
- Validated according to EN ISO 17665 / AAMI TIR 17665 (valid commissioning and product specific performance evaluation)
- Minimum sterilization temperature 132 °C (270°F); plus tolerance (0°/±3°) EN ISO 17665 / AAMI TIR 17665
- Following sterilization times (exposure time during sterilization temperature) incl. retention times are permitted:
 - Min. 5 minutes² at 134°C (273.5°F) NOT FOR US
 - Min. 3 minutes at 134°C (273.5°F) NOT FOR US
 - Min. 4 minutes at 132°C (270°F) FOR US
 - A drying time of 20 minutes is sufficient for instruments of this type

Flash sterilization is not permitted. Furthermore, do not use hot air sterilization, gamma radiation, sterilization by formaldehyde or ethylene oxide as well as plasma sterilization.

Some instruments require special aspects. Please see chapter "Special Instructions".

The fundamental suitability of an effective steam sterilization was demonstrated by an independent accredited test laboratory by application of the steam sterilizer Zentra CERT 3100 LAB (F. & M. Lautenschläger GmbH & Co. KG), Belimed Steam Sterilizer TOP 5000 Model 5-5-9, Fa. Sauter AG (K060337) The procedure described above, taking into account the worst case parameters, was taken into account.

Aufbereitungsanweisung Instrumente / Reprocessing Instructions Instruments / Instructions de préparation d'instruments / Istruzioni di purificazione degli strumenti / Instruções de reprocessamento dos instrumentos / Instrucciones de tratamiento de los instrumentos / 医疗器械准备说明 / तैयारी के निर्देश उपकरण / 器材の再処理説明書 / Инструкция по обработке инструментов / Οδηγίες προετοιμασίας εργαλείων

REF X0200 LOT 2024-04-12

7. STORING

After the sterilization, the instruments must be stored in the sterilization packaging dry and free of dust.

8. MATERIAL RESISTANCE

When selecting the detergent and disinfection agent, please see to it, that the following ingredients are not included:

- Organic, mineral and oxidizing acids (slightly alkaline detergent is recommended, pH=10 to 11, e.g. neodisher Mediclean forte from Dr. Weigert)
- Organic solvents (e.g. alcohols, ether, ketones, petrol)
- Oxidizing agents (e.g. hydrogen peroxide)
- Halogens (chlorine, iodine, bromine)
- Aromatic/halogenated hydrocarbons
- Salts of heavy metals

For cleaning all instruments, sterilization trays and sterilization containers, never use metal brushes or steel wool. All instruments, sterilization trays and sterilization containers are only allowed for temperature exposures up to 141°C (285.8°F).

For some instruments, special aspects are necessary; please see chapter "special indications"!

9. REUSABILITY

The instruments can – with appropriate care and provided they are undamaged and uncontaminated – be reused up to 200 times. Each further use exceeding this number or the use of damaged and/or contaminated instruments, is the user's responsibility. In case of non-observance, any liability is excluded.

10. RETURN

All instruments and sterilization containers must be returned reprocessed with a corresponding certificate to NGMedical or the distributor.

11. SPECIAL INSTRUCTIONS

Prior to cleaning/disinfection, all instruments must be disassembled as far as possible. In this regard, please observe the surgical technique of each system and dismantle the instruments in reverse order to the description. Impactor caps of hammers must be disassembled. Products with joints (forceps) need to be pre-cleaned, cleaned, disinfected, packed and sterilized in opened position. In case of packaging prior to sterilization, all instruments are placed into the intended holders – if available.

12. FURTHER INFORMATION

For further information, please contact

NGMedical GmbH
Trierer Straße 65
66620 Nonnweiler
GERMANY
mail@ngmedical.de
www.ngmedical.de
+49 (0) 6873 99997-0

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