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PLIF Implantate

REF X0009 LOT 2026-07-24

DE

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EN – NOT FOR US

■ INTENDED PURPOSE

BEE® PLIF cages and NGM WAVE® PLIF cages are implants for intervertebral, posterior, lumbar fusion in the segments L2-S1 among patients with a mature skeletal system.

■ PRODUCT DESCRIPTION

BEE® PLIF cages and NGM WAVE® PLIF cages have been specially adapted to suit the local anatomy to ensure the best possible surgical outcome. BEE® PLIF cages are made of titanium alloy (Ti6Al4V ELI) in accordance with ASTM F3001 and are produced in an additive manufacturing process. NGM WAVE® PLIF cages are made of PEEK-OPTIMA™ Natural and tantalum (pins) in accordance with ASTM F560. The implants are available in various heights, lengths and lordosis angles. The cranial and caudal side of the implant are rough and featurespikes to achieve a good anchoring of the implant.

■ SCOPE

These instructions for use are valid for all BEE® PLIF and NGM WAVE® PLIF implants.

■ INTENDED USERS

These instructions for use alone do not provide sufficient information for immediate use of the system. Training in the use of the products by an experienced and trained surgeon and/or by NGMedical (or the distributor) is strongly recommended.

The procedure must be performed in accordance with the instructions for use and the surgical technique. The surgeon is responsible for ensuring that the operation is performed correctly. It is strongly recommended that the procedure is only performed by surgeons who are qualified and experienced in spinal surgery, are aware of the general risks associated with spinal surgery and are familiar with the product-specific surgical techniques.

■ APPROVED INSTRUMENTS

Only the instruments described in the BEE® PLIF and NGM WAVE® PLIF surgical techniques may be used in direct connection with the implant. The surgical technique is available from NGMedical on request.

■ EXPECTED CLINICAL BENEFIT

The BEE® PLIF and NGM WAVE® lumbar cages are intended for use in lumbar fusion of adjacent vertebral bodies in accordance with their intended purpose and instructions for use. The expected clinical benefits of surgery with BEE® PLIF and NGM WAVE® include a low complication rate (especially subsidence rate) and good osseointegration (fusion) due to optimized design features such as shape and anchoring, resulting in satisfactory radiological results. Furthermore, improvement in pain and quality of life for the patient is expected. The clinical benefits are supported by preclinical data and study data.

The Summary of Safety and Clinical Performance (SSCP) can be found under the following link: www.ngmedical.de/sscp/bee_plif, www.ngmedical.de/sscp/ngm_wave

Access to the public website of the Eudamed database: <https://ec.europa.eu/tools/eudamed>.

■ INDICATIONS

BEE® PLIF and NGM WAVE® PLIF are intended for the following indications:

- Degenerative disorders of the spine including degenerative disc disease, vertebral canal and foraminal stenosis,
- Lumbar spine instability in one or more segments,
- Pseudarthrosis,
- Spondylolisthesis, including degenerative spondylolisthesis,
- Spondylodiscitis.

■ CONTRAINDICATIONS

BEE® PLIF and NGM WAVE® PLIF must not be implanted in patients with the following health conditions:

- Obesity from grade II (BMI value ≥ 35),
- Active systemic infection,
- Allergy to the implant material (see the imprint on the packaging),
- Endocrine or metabolic illness in the patient's history which influences the bone and mineral metabolism,
- Health conditions which oppose to the potential benefit of a spinal surgery,
- The patient has received medication influencing bone and mineral metabolism within two weeks prior to the planned date of operation,
- Osteoporosis,
- The patient is not willing or able to follow the post-operative instructions,
- Pregnancy,
- Vertebral fracture(s),
- Spinal metastases (exception: additional surgical stabilization measures).

■ COMPLICATIONS AND POSSIBLE UNDESIED OCCURANCES

Before the operation, the possible adverse effects and possibly required further surgery should be explained to the patient:

- Leakage of cerebrospinal fluid,
- Spinal curvature with kyphotic or scoliotic deformity, loss of correction,
- (Accelerated) adjacent segment degeneration,
- Breakage/failure of system components (loosening, degradation or breakage of components),
- Permanent disability,
- Dislocation/migration of the implant,
- Subsidence of the implant into the vertebral body end plates causing reduced height of the intervertebral space,
- Failed fusion with pseudarthrosis,
- Foreign body reactions to components, their debris or their fragments,
- Infection,
- Loss of bone density due to stress shielding,
- Possible need for further surgery,
- Nerve damage and/or spinal cord damage including the loss of neurological functions (sensory and/or motor deficiencies such as hypoesthesia, anesthesia, dysesthesia, hyperesthesia, paresthesia, partial or complete muscular paralysis, lack of reflexes, increased or reduced muscle tone, spasticity),
- Osteoporotic vertebral fracture,
- Spinal cord injuries with partial or complete paraplegia,
- Death,
- Visceral or vascular injuries,
- Spinal infections (e.g. spondylitis, discitis, myelomeningitis)
- Wound healing disorders.

■ WARNINGS

The safety and effectiveness of the BEE® PLIF and NGM WAVE® PLIF implants have only been proven for the spinal disorders listed in the "Indications" section. As with all operations, the use of this system involves the potential risk of death. Other potential risks which may require further surgery are as follows:

- Allergic reaction or sensitivity to the implant material (see the imprint on the packaging).
- Breakage/failure of system components,
- Vascular injuries,
- Pseudarthrosis,
- Damage to the spinal cord and/or nerves,
- Loss of anchorage.

The implantation of metals and alloys into the human body, exposes them to aggressive chemical environments of salts, acids and alkalis which may lead to corrosion. If different metals are positioned in close vicinity, galvanic effects may accelerate the corrosion process. In particular, BEE® PLIF and NGM WAVE® PLIF must not be used with components (exception: bone substitute material or fillers) or instruments from other manufacturers. In general, implanted material may cause immune reactions or foreign body reactions as well as chronic inflammatory changes.

When removing the implant from the packaging, it must be ensured that the nominal height of the implant matches the size stated on the packaging label.

■ PRECAUTIONS

The surgeon must be completely familiar with BEE® PLIF and NGM WAVE® PLIF, the required instruments and the general surgical technique. The surgeon must use the instruments which are specifically approved by NGMedical for use with BEE® PLIF and NGM WAVE®. The surgeon must always ensure that the adequate implant size is used and that the implant is fitted into the anatomically correct position. Overdistraction must be avoided. It must be observed that the implant is aligned correctly in the cranial-caudal as well as the ventral-dorsal plane. The general surgical procedure is described in detail in the surgical technique provided by NGMedical.

Documentation of the used implants is necessary for each patient with REF, LOT and UDI so that the level of traceability required by law is guaranteed. The implants are intended for single use. **An implant must not be reused once explanted.**

■ MRI COMPATIBILITY

The effects of magnetic resonance (imaging) on BEE® PLIF and NGM WAVE® PLIF implants have not been tested. Based on the widespread use of implants made of titanium alloy (Ti6Al4V ELI) and PEEK-OPTIMA™ Natural, BEE® PLIF and NGM WAVE® PLIF implants can be classified as conditionally MRI compatible. Patients with these implants can be scanned under 1.5 Tesla and 3.0 Tesla static magnetic field strength. A maximum specific whole-body

absorption rate (SAR) of 1.75 W/kg at 1.5 Tesla and 3.0 Tesla for a scan duration of 15 minutes should not be exceeded.

■ PACKAGING

BEE® PLIF and NGM WAVE® PLIF implants are generally supplied in double sterile packaging. Only use the implants if the packaging box is intact. Storage conditions must guarantee the integrity of the packaging box and the implants respectively. Defective implants must not be used. Each packaging has a documentation label for free use and an implant card. The packaging and labels contain all the required information. There, you can find the product number REF, LOT, description, method of sterilization, quantity, expiration date, manufacturing date and manufacturer, as well as scan the UDI code. All BEE® PLIF and NGM WAVE® PLIF implants are sterilized using gamma radiation or ethylene oxide gas. The type of sterilization is indicated on each packaging box.

■ COMPLAINTS

Every healthcare specialist (e.g. a surgeon using the product) who has a complaint or who is dissatisfied with the quality, labeling, reliability, safety, efficacy and/or performance of the product should inform NGMedical or, where applicable, the distributor. In the case of a severe incident or risk of a severe incident, liable to result in or to have resulted in the death or severe deterioration in the state of health of a patient or user or other persons, NGMedical (or the distributor) and the competent authority must be informed as soon as possible in writing or verbal. All complaints must include the product name, REF, LOT, UDI and description as well as the data of the incident. The reporting person should state his name and contact information and describe the incident in as much detail as possible. The legal reporting obligations must also be observed!

■ DISPOSAL

For disposal of NGMedical implants, please comply with the specific hospital regulations or return the processed, sterilized and labeled implants to NGMedical.


■ EXPECTED LIFETIME

The expected lifetime of BEE® PLIF and NGM WAVE® is 2.5 years after implantation under non-fused conditions. After successful fusion, the expected lifetime corresponds to the remaining lifetime of the patient. If the fusion is not successful, a revision may become necessary.

■ FURTHER INFORMATION

For further information, please contact:

 NGMedical GmbH
Trierer Straße 65
66620 Nonnweiler
GERMANY
mail@ngmedical.de
www.ngmedical.de
+49 (0) 6873 99997-0

Symbol	Standard	Ref no	Title	Meaning
	ISO 15223-1- Symbols to be used with information to be supplied by the manufacturer - Part 1: General requirements	5.4.3	Consult instructions for use	Indicates the need for the user to consult the instructions for use.
		5.1.4	Use by	Indicates the date after which the medical device is not to be used.
		5.4.2	Do not re-use	Indicates a medical device that is intended for one use, or for use on a single patient during a single procedure.
		5.2.6	Do not resterilise	Indicates a medical device that is not to be resterilised.
		5.1.6	Catalogue or model number	Indicates the manufacturer's catalogue number so that the medical device can be identified.
		5.2.8	Do not use if package is damaged	Indicates a medical device that should not be used if the package has been damaged or opened.
		5.2.4	Sterilised using irradiation	Indicates a medical device that has been sterilised using irradiation.
		5.2.3	Sterilised with ethylene oxide	Indicates a medical device that was sterilised with ethylene oxide
		5.1.1	Manufacturer	Indicates the medical device manufacturer.
		5.3.7	Storage temperature range	Indicates the temperature limits to which the medical device can be safely exposed.
		5.3.4	Keep dry	Indicates a medical device that needs to be protected from moisture.
		5.1.5	Batch code	Indicates the manufacturer's batch code so that the batch or lot can be identified.
		5.1.3	Date of manufacture	Indicates the date of manufacture the medical device
		5.7.7	Medical Device	Indicates the item is a medical device
		5.7.10	Unique Device Identification	Indicates a carrier that contains Unique Device Identifier information
	Double sterile barrier system			
	European conformity mark with Notified Body identification number			
	The medical device according to 21 CFR 801.109(b)(1) may only be prescribed by a physician in accordance with the law (USA)			
	Conditionally MRI compatible			

® The NGMedical brands and registered trademarks are protected by national and international rights in Germany, the European Union and other regions and countries around the world.

INSTRUCTIONS FOR USE

PLIF implants

REF **X0009** **LOT** 2026-07-24

EN-FOR US ONLY

■ CAUTION

Federal law (USA) restricts these devices to sale by or on the order of a physician. These devices should be implanted only by a physician who is fully trained with the devices, intended use, instrumentation and with knowledge of the surgical techniques required. Contact your representative for surgical technique.

■ INTENDED PURPOSE

BEE® PLIF cages and NGM WAVE® PLIF cages are implants for intervertebral, posterior, lumbar fusion in the segments L2-S1 among patients with a mature skeletal system.

■ PRODUCT DESCRIPTION

BEE® PLIF cages and NGM WAVE® PLIF cages have been specially adapted to suit the local anatomy to ensure the operative outcome to the greatest extent possible. BEE® PLIF cages are made of titanium alloy in accordance with ASTM F3001 and are produced in an additive manufacturing process. NGM WAVE® PLIF cages are made of PEEK-OPTIMA™ Natural and tantalum (pins) in accordance with ASTM F560. The implants are available in various heights, lengths and lordosis angles. The cranial and caudal side of the implant are rough and feature spikes to achieve a good degree of anchoring of the implant.

■ SCOPE

These instructions for use are valid for all BEE® PLIF implants and NGM WAVE® PLIF.

■ APPROVED INSTRUMENTS

Only the instruments described in the BEE® PLIF surgical techniques may be used in direct connection with the implant. The removal of the implant is described in the surgical technique, which you can obtain from NGMedical upon request.

■ INDICATIONS

BEE® PLIF cages and NGM WAVE® PLIF cages are indicated for interbody fusion with autogenous bone graft in patients with degenerative disc disease (DDD) at one or two contiguous levels from L2 to S1. These DDD patients may also have up to Grade 1 spondylolisthesis or retrolisthesis at the involved levels. DDD is defined as discogenic back pain with degeneration of the disc confirmed by history and radiographic studies. These patients should be skeletally mature and have had six months of non-operative treatment. These implants may be implanted via an open or a minimally invasive posterior approach. These implants are to be used with autogenous and/or allogeneic bone graft comprised of cancellous and/or corticocancellous bone graft to facilitate fusion. These devices are intended to be used with supplemental fixation instrumentation, which has been cleared by the FDA for use in the lumbar spine.

■ CONTRAINDICATIONS

BEE® PLIF and NGM WAVE® PLIF cages must not be implanted in patients with the following health situations:

- Obesity from grade II (BMI value ≥ 35),
- Active system infection,
- Allergy to the implant material (see the imprint on the packaging).
- Endocrine or metabolic illness in the patient's history which influences the bone and mineral balance,
- Health conditions which oppose the potential benefit of a spinal surgery,
- The patient took medication which influences the bone and mineral balance within two weeks prior to the planned date of operation,
- Osteoporosis,
- The patient is not willing or able to follow the post-operative instructions,
- Pregnancy,
- Vertebral fracture(s),
- Spinal metastases (exception: additional surgical stabilization measures).

■ COMPLICATIONS AND POSSIBLE UNDESIRABLE OCCURRENCES

Before the operation, the possible adverse effects and further surgery which may be necessary to correct the effects should be explained to the patient:

- Discharge of cerebrospinal fluid,
- Change to the spinal curvature with kyphotic or scoliotic deformity, loss of correction,
- (Accelerated) adjacent segment degeneration,
- Breakage/failure of system components (loosening, degradation or breakage of components),
- Permanent disability,
- Dislocation/migration of the implant,
- Subsidence of the implant into the vertebral body end plates with a reduced height of the intervertebral space,
- Failed fusion with pseudarthrosis,
- Foreign body granuloma around components, their debris or their fragments,
- Infection,
- Loss of bone density due to stress shielding,
- Possible need for further surgery,
- Nerve damage and/or spinal cord damage including the loss of neurological functions (sensory and/or motor deficiencies such as hyposthesia, anesthesia, dyesthesia, hyperesthesia, paresthesia, partial or complete muscular paralysis, lack of reflexes, increased or reduced muscle tone, spasticity),
- Osteoporotic vertebral fracture,
- Spinal cord injuries with partial or complete paraplegia,
- Death,
- Visceral or vascular injuries,
- Spinal infections (e.g. spondylitis, discitis, myelomeningitis)
- Wound healing disorders.

■ WARNINGS

The safety and effectiveness of the BEE® PLIF and NGM WAVE® PLIF implants have only been proven for the spinal disorders listed in the "Indications" section. As with all operations, the use of this system involves the potential risk of death. Other potential risks which may require further surgery are as follows:

- Allergic reaction or sensitivity to the implant material (see the imprint on the packaging).
- Breakage/failure of system components,
- Vascular injuries,
- Pseudarthrosis,
- Damage to the spinal cord and/or nerves,
- Loss of anchorage.

The implantation of metals and alloys into the human body, expose them to aggressive chemical environments of salts, acids and bases which may lead to corrosion. If different metals are positioned next to one another in touching proximity, galvanic corrosion may accelerate the corrosion process. In particular, BEE® PLIF and NGM WAVE® PLIF must not be used with components (exception: autogenous and/or allogeneic bone graft comprised of cancellous and/or corticocancellous bone graft to facilitate fusion) or instruments from other manufacturers. In general, implanted material may cause immune reactions or foreign body reactions as well as chronic inflammatory changes.

When removing the implant from the packaging, it must be ensured that the nominal height of the implant matches the size stated on the packaging label.

■ PRECAUTIONS

The surgeon must be completely familiar with BEE® PLIF and NGM WAVE® PLIF, the required instruments and the general surgical technique. The surgeon must use the instruments which are specifically approved by NGMedical for use with BEE® PLIF and NGM WAVE® PLIF. Based on the dynamic testing results, the physician should consider the levels of implantation, patient weight, patient activity level, other patient conditions, etc., which may impact on the performance of the intervertebral body fusion device. The surgeon must always ensure that the ideal implant size is used and that the implant is fitted into the anatomically correct position. Overdistraction must be avoided. Here, it must be observed that the implant is aligned correctly both in terms of the cranial-caudal line and the ventral-dorsal line. The general surgical procedure is described in detail in the surgical technique provided by NGMedical.

Please document the used implants for each patient with REF, LOT and UDI so that the level of traceability required by law is guaranteed. The implants are intended for single use. **An implant must not be reused once explanted.**

■ MRI COMPATIBILITY

BEE® PLIF and NGM WAVE® PLIF have not been evaluated for safety and compatibility in the MR environment. It has not been tested for heating migration, or image artifact in the MR environment. The safety of the BEE® PLIF and NGM WAVE® PLIF systems in the MR environment is unknown. Scanning a patient who has this device may result in patient injury.

■ PACKAGING

BEE® PLIF and NGM WAVE® PLIF implants are generally supplied in double sterile packaging. Only use the implants if the packaging is intact. The storage conditions must guarantee the intact nature of the packaging and therefore also the implants. Defective implants must not be used. Each packaging has a documentation label for you to use and an implantation ID. The packaging and labels contain all the required information. There, you can find the product number REF, LOT, description, method of sterilization, quantity, expiration date, manufacturing date and manufacturer, as well as scan the UDI code. All BEE® PLIF and NGM WAVE® PLIF implants are sterilized by way of gamma irradiation or EO gas sterilization. The type of sterilization is indicated on each packaging box.

■ COMPLAINTS

Every healthcare specialist (e.g. a surgeon using the product) who has a complaint or who is dissatisfied with the quality, labeling, reliability, safety, efficacy and/or performance of the product should inform NGMedical or, where

applicable, the distributor. In the case of a severe incident or the risk of a severe incident, liable to result in or to have resulted in the death or severe deterioration in the state of health of a patient or user or other persons, NGMedical (or the distributor) and the competent authority must be informed as soon as possible in writing or verbal. All complaints must include the product name, REF, LOT, UDI and the description as well as the data of the incident. The reporting person should state his name, contact data and a detailed description of the incident. The legal reporting obligations are also to be observed!

■ DISPOSAL

For disposal of NGMedical implants, please comply with the specific hospital regulations or return the processed, sterilized and labeled implants to NGMedical.


■ EXPECTED LIFETIME

The expected lifetime of BEE® PLIF and NGM WAVE® is 2.5 years after implantation under non-fused conditions. After successful fusion, the expected lifetime corresponds to the remaining lifetime of the patient. If the fusion is not successful, a revision may become necessary.

■ FURTHER INFORMATION

For further information, please contact:

 NGMedical GmbH
Trierer Straße 65
66620 Nonnweiler
GERMANY
mail@ngmedical.de
www.ngmedical.de
+49 (0) 6873 99997-0

The Summary of Safety and Clinical Performance (SSCP) can be found under the following link:
www.ngmedical.de/sscp/bee_plif, www.ngmedical.de/sscp/ngm_wave

Access to the public website of the Eudamed database: <https://ec.europa.eu/tools/eudamed>.

Symbol	Standard	Ref no	Title	Meaning
	ISO 15223-1- Symbols to be used with information to be supplied by the manufacturer - Part 1: General requirements	5.4.3	Consult instructions for use	Indicates the need for the user to consult the instructions for use.
		5.1.4	Use by	Indicates the date after which the medical device is not to be used.
		5.4.2	Do not re-use	Indicates a medical device that is intended for one use, or for use on a single patient during a single procedure.
		5.2.6	Do not re-sterilize	Indicates a medical device that is not to be re-sterilized.
		5.1.6	Catalogue or model number	Indicates the manufacturer's catalogue number so that the medical device can be identified.
		5.2.8	Do not use if package is damaged	Indicates a medical device that should not be used if the package has been damaged or opened.
		5.2.4	Sterilized using irradiation	Indicates a medical device that has been sterilized using irradiation.
		5.2.3	Sterilized using EO	Indicates a medical device that has been sterilized using EO gas.
		5.1.1	Manufacturer	Indicates the medical device manufacturer.
		5.3.7	Storage temperature range	Indicates the temperature limits to which the medical device can be safely exposed.
		5.3.4	Keep dry	Indicates a medical device that needs to be protected from moisture.
		5.1.5	Batch code	Indicates the manufacturer's batch code so that the batch or lot can be identified.
		5.1.3	Date of manufacture	Indicates the date of manufacture the medical device
		5.7.7	Medical Device	Indicates the item is a medical device
		5.7.10	Unique Device Identification	Indicates a carrier that contains Unique Device Identifier information
	Double sterile barrier system			
	European conformity mark with Notified Body identification number			
	The medical device according to 21 CFR 801.109(b)(1) may only be prescribed by a physician in accordance with the law (USA)			
	Conditionally MRI compatible			

© The NGMedical brands and registered trademarks are protected by national and international rights in Germany, the European Union and other regions and countries around the world.

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