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EN | NOT FOR US

■ INTENDED PURPOSE

The ART® spinal fixation system is intended for the posterior stabilization of the thoracic, lumbar and lumbosacral spine among patients with a mature skeletal system.

■ PRODUCT DESCRIPTION

The ART® spinal fixation system was developed with the goal of providing the surgeon with a fixation system to carry out dorsal spinal stabilizations simply, quickly and effectively.

Important results of these developments are:

- Polyaxial pedicle screw with rod connector (screw head)
- Pedicle screw with double thread
- Freely rotatable rod connector
- Top-loading system
- Option of adjusting the rod connector by up to ± 30° (60°)

ART® is made from a titanium alloy (Ti6Al4V_ELI) in accordance with the standard ASTM F136. ART® must not be used in direct connection with components from other manufacturers.

■ SCOPE

These instructions for use are valid for all implants of the ART® product group.

■ INTENDED USERS

These instructions for use alone do not provide sufficient information for the immediate use of the system. Training in the use of the products by an experienced and trained surgeon and/or by NGMedical (or the distributor) is strongly recommended.

The procedure must be performed in accordance with the instructions for use and the surgical technique. The surgeon is responsible for ensuring that the operation is performed correctly. It is strongly recommended that the procedure is only performed by surgeons who are qualified and experienced in spinal surgery, are aware of the general risks associated with spinal surgery and are familiar with the product-specific surgical techniques.

■ APPROVED INSTRUMENTS

Only the instruments described in the ART® surgical technique may be used in direct connection with the implant. The surgical technique is available from NGMedical on request.

■ EXPECTED CLINICAL BENEFIT

The ART® spinal fixation system provides stabilisation of the spine in accordance with its intended purpose and instructions for use. The expected clinical benefit of a surgery with ART® include low complication rates, quick and easy implantation, and satisfactory radiological results. This is expected to result in pain reduction and support fusion for the patient.

The Summary of Safety and Clinical Performance (SSCP) report is available at the following link: www.ngmedical.de/sscp/art. Access to the public website of the Eudamed database: <https://ec.europa.eu/tools/eudamed>.

■ INDICATIONS

The ART® spinal fixation system is intended to be used for the following indications:

- To fix and stabilize vertebral segments,
 - To support osseous fusion,
 - When treating the following acute and chronic instabilities and deformities of the thoracic, lumbar and sacral spine: degenerative spondylolisthesis with objective evidence of a relevant spinal stenosis, fractures, dislocations, scoliosis, kyphosis, spinal tumors and failed prior fusions (pseudoarthrosis),
 - When treating severe spondylolisthesis (grade 3 and 4) of the lumbar and lumbosacral vertebral segments among patients with a mature skeletal system.
- Furthermore, ART® can also be used in connection with MOVE®-P. Please refer to the MOVE®-P instructions for use and surgical technique.

■ CONTRAINDICATIONS

- Obesity from grade II (BMI value ≥ 35),
 - Active systemic infection,
 - Allergy to the implant material (Ti6Al4V_ELI),
 - Endocrine or metabolic illnesses in the patient's history which influence the bone and mineral metabolism,
 - Health conditions which are opposed to the potential benefit of a spinal operation,
 - The patient has received medication influencing bone and mineral metabolism less than two weeks before the planned date of operation,
 - Bone anomalies which prevent the secure fixation of screws (e.g. dysplastic pedicles),
 - The patient is not willing or able to follow the post-operative instructions,
 - Pregnancy,
 - Spinal metastases (exception: additional surgical stabilization measures).
- ART® was neither developed for indications other than those stated above, nor is it marketed for these.

■ COMPLICATIONS AND POSSIBLE UNDESIRE OCCURENCES

Before the operation, possible adverse effects and further surgery which may be necessary to correct these effects should be explained to the patient:

- Adjacent segment degeneration,
- Leakage of cerebrospinal fluid,
- Failed fusion,
- Foreign body reactions to components (caused by allergies), their debris or their fragments,
- Infections, including those of the meninges,
- Loosening, detachment or breakage of components,
- Nerve damage and/or spinal cord damage including the loss of neurological functions (sensory and/or motor deficits such as hypoesthesia, anesthesia, dysesthesia, hyperesthesia, paresthesia, muscular paresis or paralysis, loss of reflexes, increased or reduced muscular tonus, spasticity, bladder, bowel and sexual dysfunctions),
- Necessity of further surgery,
- Spinal cord injuries with partial or complete paraplegia,
- Death,
- Visceral or vascular injuries,
- Wound healing disorders,
- Osteolysis.

■ WARNINGS

The safety and effectiveness of the ART® implants has been proven only for the spinal disorders stated in the indications section. As with all operations, the use of this system involves the potential risk of death. Other potential risks which may require additional surgery are as follows:

- Allergic reaction or sensitivity to the titanium alloy (Ti6Al4V_ELI),
- Deformation and/or breakage of system components,
- Infections,
- Loosening of components,
- Pseudoarthrosis,
- Damage to the spinal cord and/or nerves.

By implantation of metals and alloys into the human body, these are exposed to aggressive chemical milieus of salts, acids and alkalis, which might lead to corrosion. If different metals are positioned in close vicinity, galvanic effects may accelerate the corrosion process. In particular, ART® must not be used with components (exception: bone substitute material or fillers) or instruments from other manufacturers.

When removing the implant from the packaging, it must be ensured that the nominal size of the implant matches the size stated on the packaging label.

A completely satisfactory outcome cannot be achieved in every operation. This also applies to spinal operations, where the outcome is determined by many factors.

The osteosynthesis with ART® must normally be supported by an arthrodesis of the affected vertebrae.

■ PRECAUTIONS

The surgeon must be perfectly familiar with ART®, the required instruments and the general surgical technique. The surgeon must exclusively use the instruments especially approved by NGMedical for ART®. For maximum stability, the surgeon must select the correct size for the pedicle. The general surgical procedure is described in detail in the surgical technique provided by NGMedical (REF X1006).

Please document the used implants for each patient with REF, LOT and UDI so that the level of traceability required by law is guaranteed. The implants are intended for single use. **An explanted implant must never be reused.**

■ MRI COMPATIBILITY

The effects of magnetic resonance (imaging) on ART® implants have not been tested. Based on the widespread use of implants made of titan alloy (Ti6Al4V_ELI), ART® implants can be classified as conditionally MRI compatible. Patients with these implants can be scanned under 1.5 Tesla and 3.0 Tesla static magnetic field strength. A maximum specific whole-body absorption rate (SAR) at 1.5 Tesla and 3.0 Tesla of 1.75 W/kg with a scan duration of 15 minutes should not be exceeded.

■ PACKAGING

ART® implants are generally supplied in double sterile packaging. Please only use implants if the packaging box is undamaged. The storage conditions must guarantee the integrity of the packaging box and the implants respectively. Defective implants or other components must not be used. Each packaging box contains documentation labels for free use and an implant card. The packaging box and labels contain all information required. They indicate the product number REF, LOT, description, sterilization method, quantity, expiration date, manufacturing date and manufacturer as well as the UDI code. All ART® implants are sterilized using gamma radiation or ethylene oxide gas. The type of sterilization is indicated on each packaging box.

■ COMPLAINTS

Any healthcare professional (e.g. surgeon using the product) who has a complaint or who is dissatisfied with the quality, identification, reliability, safety, efficacy and/or performance of the ART® products, should inform NGMedical or, where applicable, the distributor. In the case of a severe incident or the risk of a severe incident, liable to result in or to have resulted in the death or severe deterioration in the state of health of a patient or a user or other persons, NGMedical (or the distributor) and the competent authority must be informed as soon as possible in writing or verbal. All complaints must include the product name, REF, LOT, UDI and description as well as the data of the incident. The reporting person should state his name, contact data and a detailed description of the incident. The legal reporting obligations are also to be observed!

■ DISPOSAL

For disposal of NGMedical implants, please comply with the specific hospital regulations or return the processed, sterilized and labeled implants to NGMedical.


■ EXPECTED LIFETIME

The expected lifetime of ART® is 2 years after implantation without fusion. After successful fusion, the expected lifetime corresponds to the remaining lifetime of the patient. In combination with the motion-preserving implant MOVE®-P by NGMedical GmbH, the expected lifetime of the combined system is adapted to the shortest lifetime of the individual implants used. As the expected lifetime of MOVE®-P is designed for 10 years, the expected lifetime of the combination of ART® and MOVE®-P is also 10 years. After exceeding the expected lifetime, a revision may become necessary.

■ FURTHER INFORMATION

For further information, please contact:

NGMedical GmbH
Trierer Straße 65
66620 Nonnweiler
GERMANY
mail@ngmedical.de
www.ngmedical.de
+49 (0) 6873 99997-0

Symbol	Standard	Ref no	Title	Meaning
		5.4.3	Consult instructions for use	Indicates the need for the user to consult the instructions for use.
		5.1.4	Use by	Indicates the date after which the medical device is not to be used.
		5.4.2	Do not re-use	Indicates a medical device that is intended for one use, or for use on a single patient during a single procedure.
		5.2.6	Do not re-sterilize	Indicates a medical device that is not to be re-sterilized.
		5.1.6	Catalogue or model number	Indicates the manufacturer's catalogue number so that the medical device can be identified.
		5.2.8	Do not use if package is damaged	Indicates a medical device that should not be used if the package has been damaged or opened.
	ISO 15223-1- Symbols to be used with information to be supplied by the manufacturer - Part 1: General requirements	5.2.4	Sterilized using irradiation	Indicates a medical device that has been sterilized using irradiation.
		5.2.3	Sterilized with ethylene oxide	Indicates a medical device that was sterilized with ethylene oxide
		5.1.1	Manufacturer	Indicates the medical device manufacturer.
		5.3.7	Storage temperature range	Indicates the temperature limits to which the medical device can be safely exposed.
		5.3.4	Keep dry	Indicates a medical device that needs to be protected from moisture.
		5.1.5	Batch code	Indicates the manufacturer's batch code so that the batch or lot can be identified.
		5.1.3	Date of manufacture	Indicates the date of manufacture the medical device
		5.7.7	Medical Device	Indicates the item is a medical device
		5.7.10	Unique Device Identification	Indicates a carrier that contains Unique Device Identifier information
	Double sterile barrier system			
	European conformity mark with Notified Body identification number			
	The medical device according to 21 CFR 801.109(b)(1) may only be prescribed by a physician in accordance with the law (USA)			
	Conditionally MRI compatible			

© The NGMedical brands and registered trademarks are protected by national and international rights in Germany, the European Union and other regions and countries around the world.

INSTRUCTIONS FOR USE

ART® Implants

REF X0006 **LOT** 2025-08-29

EN | FOR US ONLY

■ CAUTION

The Federal law (USA) restricts these devices to sale by or on the order of a physician. These devices should be implanted only by a physician who is fully trained with the devices, intended use, instrumentation and with knowledge of the surgical techniques required. Contact your representative for surgical technique.

■ INTENDED PURPOSE

The ART® spinal fixation system is intended for the posterior stabilization of the thoracic, lumbar and lumbosacral spine among patients with a mature skeletal system.

■ PRODUCT DESCRIPTION

The ART® spinal fixation system was developed with the goal of providing the surgeon with a fixation system to carry out dorsal spinal stabilizations simply, quickly and effectively.

Important results of these developments are:

- Polyaxial pedicle screw with rod connector (screw head)
- Pedicle screw with double thread
- Freely rotatable rod connector
- Top-loading system
- Option of adjusting the rod connector by up to ± 30° (60°)

ART® is made from a titanium alloy (Ti6Al4V_ELI) in accordance with the standard ASTM F136. ART® must not be used in direct connection with components from other manufacturers. This does not include the use of bone cement.

■ SCOPE

These instructions for use are valid for all implants of the ART® product group.

■ APPROVED INSTRUMENTS

Only the instruments described in the ART® surgical technique may be used in direct connection with the implant. The surgical technique is available from NGMedical on request.

■ INDICATIONS

The ART® spinal fixation system is intended to be used for the following indications:

- To fix and stabilize vertebral segments,
- To support osseous fusion,
- When treating the following acute and chronic instabilities and deformities of the thoracic, lumbar and sacral spine: degenerative spondylolisthesis with objective evidence of a relevant spinal stenosis, fractures, dislocations, scoliosis, kyphosis, spinal tumors and failed prior fusions (pseudoarthrosis),
- When treating severe spondylolisthesis (grade 3 and 4) of the lumbar and lumbosacral vertebral segments among patients with a mature skeletal system.

■ CONTRAINDICATIONS

- Obesity from grade II (BMI value ≥ 35),
- Active systemic infection,
- Allergy to the implant material (Ti6Al4V_ELI),
- Endocrine or metabolic illnesses in the patient's history which influence the bone and mineral metabolism,
- Health conditions which are opposed to the potential benefit of a spinal operation,
- The patient has received medication influencing bone and mineral metabolism less than two weeks before the planned date of operation,
- Bone anomalies which prevent the secure fixation of screws (e.g. dysplastic pedicles),
- The patient is not willing or able to follow the post-operative instructions,
- Pregnancy,
- Spinal metastases (exception: additional surgical stabilization measures).

ART® was neither developed for indications other than those stated above, nor is it marketed for these.

■ COMPLICATIONS AND POSSIBLE ADVERSE EFFECTS

Before the operation, possible adverse effects and further surgery which may be necessary to correct these effects should be explained to the patient:

- Adjacent segment degeneration,
- Leakage of cerebrospinal fluid,
- Failed fusion,
- Foreign body reactions to components (caused by allergies), their debris or their fragments,
- Infections, including those of the meninges,
- Loosening, detachment or breakage of components,
- Nerve damage and/or spinal cord damage including the loss of neurological functions (sensory and/or motor deficits such as hypoesthesia, anesthesia, dysesthesia, hyperesthesia, paresthesia, muscular paresis or paralysis, loss of reflexes, increased or reduced muscular tonus, spasticity, bladder, bowel and sexual dysfunctions),
- Necessity of further surgery,
- Spinal cord injuries with partial or complete paraplegia,
- Death,
- visceral or vascular injuries,
- wound healing disorders,
- osteolysis.

■ WARNINGS

The safety and effectiveness of the ART® implants has been proven only for the spinal disorders stated in the indications section. As with all operations, the use of this system involves the potential risk of death. Other potential risks which may require additional surgery are as follows:

- Allergic reaction or sensitivity to the titanium alloy (Ti6Al4V_ELI),
- Deformation and/or breakage of system components,
- Infections,
- Loosening of components,
- Pseudoarthrosis,
- Damage to the spinal cord and/or nerves.

By implantation of metals and alloys into the human body, these are exposed to aggressive chemical milieu of salts, acids and alkalis, which might lead to corrosion. If different metals are positioned in close vicinity, galvanic effects may accelerate the corrosion process. In particular, ART® must not be used with components (exception: bone substitute material or fillers) or instruments from other manufacturers.

When removing the implant from the packaging, it must be ensured that the nominal size of the implant matches the size stated on the packaging label.

A completely satisfactory outcome cannot be achieved in every operation. This also applies to spinal operations, where the outcome is determined by many factors.

The osteosynthesis with ART® must normally be supported by an arthrodesis of the affected vertebrae.

■ PRECAUTIONS

The surgeon must be perfectly familiar with ART®, the required instruments and the general surgical technique. The surgeon must exclusively use the instruments especially approved by NGMedical for ART®. For maximum stability, the surgeon must select the correct size for the pedicle. The general surgical procedure is described in detail in the surgical technique provided by NGMedical (REF X1006).

Please document the used implants for each patient with REF, LOT and UDI so that the level of traceability required by law is guaranteed. The implants are intended for single use. **An explanted implant must never be reused.**

■ MRI COMPATIBILITY

ART® implants have not been evaluated for safety and compatibility in the MR environment. They have not been tested for heating migration, or image artifact in the MR environment. The safety of the ART® system in the MR environment is unknown. Scanning a patient who has this device may result in patient injury.

■ PACKAGING

ART® implants are generally supplied in double sterile packaging. Please only use implants if the packaging box is undamaged. The storage conditions must guarantee the integrity of the packaging box and the implants respectively. Defective implants or other components must not be used. Each packaging box contains documentation labels for free use and an implant card. The packaging box and labels contain all information required. They indicate the product number REF, LOT, description, sterilization method, quantity, expiration date, manufacturing date and manufacturer as well as the UDI code. All ART® implants are sterilized using gamma radiation or ethylene oxide gas. The type of sterilization is indicated on each packaging box.

■ COMPLAINTS

Any healthcare professional (e.g. surgeon using the product) who has a complaint or who is dissatisfied with the quality, identification, reliability, safety, efficacy and/or performance of the ART® products, should inform NGMedical or, where applicable, the distributor. In the case of a severe incident or the risk of a severe incident, liable to result in or to have resulted in the death or severe deterioration in the state of health of a patient or a user or other persons, NGMedical (or the distributor) and the competent authority must be informed as soon as possible in writing or verbal. All complaints must include the product name, REF, LOT, UDI and description as well as the date of the incident. The reporting person should state his name, contact data and a detailed description of the incident. The legal reporting obligations are also to be observed!

■ DISPOSAL

For disposal of NGMedical implants, please comply with the specific hospital regulations or return the processed, sterilized and labeled implants to NGMedical.


■ EXPECTED LIFETIME

The expected lifetime of ART® is 2 years after implantation without fusion. After successful fusion, the expected lifetime corresponds to the remaining lifetime of the patient. In direct combination with MOVE-P, the expected lifetime is 10 years. After exceeding the expected lifetime, a revision may become necessary.

■ FURTHER INFORMATION






For further information, please contact:

 NGMedical GmbH
Trierer Straße 65
66620 Nonnweiler
GERMANY
mail@ngmedical.de
www.ngmedical.de
+49 (0) 6873 99997-0

The Summary of Safety and Clinical Performance (SSCP) can be found under the following link: www.ngmedical.de/sscp/art.

Access to the public website of the Eudamed database: <https://ec.europa.eu/tools/eudamed>.

Symbol	Standard	Ref no	Title	Meaning
		5.4.3	Consult instructions for use	Indicates the need for the user to consult the instructions for use.
		5.1.4	Use by	Indicates the date after which the medical device is not to be used.
		5.4.2	Do not re-use	Indicates a medical device that is intended for one use, or for use on a single patient during a single procedure.
		5.2.6	Do not re-sterilize	Indicates a medical device that is not to be re-sterilized.
		5.1.6	Catalogue or model number	Indicates the manufacturer's catalogue number so that the medical device can be identified.
		5.2.8	Do not use if package is damaged	Indicates a medical device that should not be used if the package has been damaged or opened.
	ISO 15223-1- Symbols to be used with information to be supplied by the manufacturer - Part 1: General requirements	5.2.4	Sterilized using irradiation	Indicates a medical device that has been sterilized using irradiation.
		5.2.3	Sterilized with ethylene oxide	Indicates a medical device that was sterilized with ethylene oxide.
		5.1.1	Manufacturer	Indicates the medical device manufacturer.
		5.3.7	Storage temperature range	Indicates the temperature limits to which the medical device can be safely exposed.
		5.3.4	Keep dry	Indicates a medical device that needs to be protected from moisture.
		5.1.5	Batch code	Indicates the manufacturer's batch code so that the batch or lot can be identified.
		5.1.3	Date of manufacture	Indicates the date of manufacture the medical device
		5.7.7	Medical Device	Indicates the item is a medical device
		5.7.10	Unique Device Identification	Indicates a carrier that contains Unique Device Identifier information
	Double sterile barrier system			
	European conformity mark with Notified Body identification number			
	The medical device according to 21 CFR 801.109(b)(1) may only be prescribed by a physician in accordance with the law (USA)			
	Conditionally MRI compatible			

© The NGMedical brands and registered trademarks are protected by national and international rights in Germany, the European Union and other regions and countries around the world.

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